# 

# **Spectrophotometers**





- Visible
- UV/Vis
- SpectroQuest<sup>™</sup>
- Software
- Accessories
- Cuvets







#### **Table of Contents**



elcome to the UNICO Spectrophotometer catalog. Inside you will find our complete line of spectrophotometers and accessories. Our line ranges from very basic Visible wavelength instruments, to advanced Scanning UV/Visible spectrophotometers. As you browse the catalog, we are confident you will find an instrument to fit your applications and budget.

#### **Contact UNICO:**

UNICO is headquartered in Dayton, NJ, U.S.A., just outside of Princeton, where our office, warehouse and service center are located. UNICO is open for business Monday- Friday, 52 weeks a year from 8:30 – 5:30 pm Eastern Standard Time. We are closed for certain national holidays. We can be reached by phone in the U.S. toll free at 800-588-9776; if dialing internationally use 732-274-1155. Our fax number is 732-274-1151.

Please visit our website at www.unicosci.com or e-mail us anytime at info@unicosci.com or Sales@unicosci.com.

#### **UNICO**

182 Ridge Road, Suite E Dayton, New Jersey 08810 U.S.A. www.unicosci.com

#### Ordering Information:

UNICO products are available from local, national and international distributors worldwide. If you are interested in becoming a UNICO distributor, please contact us at: info@unicosci.com for a distributorship application. If you would like to purchase UNICO products and are unable to find a local supplier, please contact UNICO at info@unicsoci.com and we will direct you to the nearest supplier.

Copyright United Products and Instruments, Inc. dba UNICO © Copyright UNICO 2009, all rights reserved

#### Estimado cliente,

Muchas gracias por su interés en nuestros productos. Este nuevo catálogo incluye todos nuestros espectrofotómetros. UNICO también ofrece microscopios, centrifugas, mezcladores, agitadores y otros equipos y accesorios de laboratorio. Para servicio en español, diríjase directamente a nuestra página Web en donde encontrará folletos de todos nuestros productos, traducidos a nuestro idioma y publicados en formato PDF.

También nos encontramos a vuestra disposición por teléfono al 609-240-5507 o por email a ventas@unicosci.com. Será un placer servirlo.

Table of Contents	Pages
Visible Spectrophotometer Selection Guide	3
Basic Visible Spectrophotometers, Series 1000	4
Visible Spectrophotometers, Series 1100	5-6
Visible Spectrophotometers, Series 1200 - 1205	7-9
Visible Spectrophotometers, Long Path, Series 2100+	10
UV/Vis Spectrophotometer Selection Guide	11
UV/Visible Spectrophotometers, Series 2100UV+	12-14
UV/Vis SpectroQuest™ Scanning Spectrophotometers, Series 2800	15-17
UV/Vis SpectroQuest™ Scanning Spectrophotometers, Series 3800	18
UV/Vis SpectroQuest™ Scanning Spectrophotometers, Series 4800	19
SpectroQuest™ Spectrophotometer Software	20-23
SpectroQuest™ Accessories	24-25
Spectrophotometer Cuvet Selection Guide	26
UNICO Laboratory Equipment	27



UNICO manufactures a complete line of Centrifuges, Microscopes, Lab Equipment, and Accessories not shown in this catalog. Please visit our website at www.unicosci.com to see our entire range of products. Downloadable PDF brochures are available on the site.

Model MX PowerSpin™ Centrifuge



Model MTR22 Programmable Multi-Mix Rotator



#### **Visible Spectro Selection Guide**

#### **Visible Spectrophotometer Selection Guide**



\$1000

#### **Electromagnetic Spectrum:**

Wavelength Range in Nanometers (nm):

Ultra-Violet (UV), Visible (Vis) and Infrared Light (IR)

Ultra-Violet (UV) Spectrum: 200nm - 400nm
 Visible Light Spectrum (Vis): 400nm - 750nm

• Infrared Spectrum (IR): 750nm – 300,000nm



S1200

#### **UNICO Visible Light Spectrophotometers**

Wavelength Range by series:

\$1000 Series: 400nm - 1,000nm
 \$1100 Series: 335nm - 1,000nm
 \$1200 Series: 325nm - 1,000nm
 \$2100+ Series: 325nm - 1,000nm



S1100 S2100+

How do I select the right spectrophotometer for my applications? We can help. Contact UNICO by phone, fax, or e-mail and our technicians will be glad to help select the right instrument for your needs. Please see the Selection Guide below. This will help you understand the general specifications by model.

Model	Wavelength Range (nm)	Slit Width Bandpass	Cell Pathlength	Scanning	Absorbance Range	Concentration Conc., C	Photometric Accuracy	Data Port(s)	Light Beam	Optional Printer	Voltage Pre-set	Pg.
<b>S1000</b>	400-1000 nm	20 nm	10 mm	N/A	0-1.99A	N/A	+/- 2%	USB	Single	N/A	110 Volt	4
S1000E	400-1000 nm	20 nm	10 mm	N/A	0-1.99A	N/A	+/- 2%	USB	Single	N/A	220 Volt	4
S1100	335-1000 nm	20 nm	10 mm	N/A	0-2.0A	N/A	+/- 2%	USB	Single	N/A	110 Volt	5
S1100E	335-1000 nm	20 nm	10 mm	N/A	0-2.0A	N/A	+/- 2%	USB	Single	N/A	220 Volt	5
S1100RS	335-1000 nm	10 nm	10 mm	N/A	0-2.0A	0 - 1999 C	+/- 1%	USB & RS232	Single	N/A	110 Volt	6
S1100RSE	335-1000 nm	10 nm	10 mm	N/A	0-2.0A	0 - 1999 C	+/- 1%	USB & RS232	Single	N/A	220 Volt	6
<b>S1200</b>	325-1000 nm	5 nm	50 mm	N/A	0-2.0A	0 - 1999 C	+/- 1%	USB & RS232	Single	N/A	110 Volt	7
S1200E	325-1000 nm	5 nm	50 mm	N/A	0-2.0A	0 - 1999 C	+/- 1%	USB & RS232	Single	N/A	220 Volt	7
S1201	325-1000 nm	5 nm	50 mm	N/A	0-2.0A	0 - 1999 C	+/- 1%	USB & RS232	Single	N/A	110 Volt	7
S1201E	325-1000 nm	5 nm	50 mm	N/A	0-2.0A	0 - 1999 C	+/- 1%	USB & RS232	Single	N/A	220 Volt	7
\$1205	325-1000 nm	5 nm	50 mm	Via PC	0.3-2.5A	0 - 9999 C	.004 @ 0.5 A	USB & RS232	Single	N/A	110 Volt	8
S1205E	325-1000 nm	5 nm	50 mm	Via PC	0.3-2.5A	0 - 9999 C	.004 @ 0.5 A	USB & RS232	Single	N/A	220 Volt	8
S2100+	325-1000 nm	4 nm	100 mm	Via PC	0.3-2.5A	0 - 9999 C	.004 @ 0.5 A	USB & RS232	Single	Yes	110 Volt	10
S2100+E	325-1000 nm	4 nm	100 mm	Via PC	0.3-2.5A	0 - 9999 C	.004 @ 0.5 A	USB & RS232	Single	Yes	220 Volt	10
S2100+P	325-1000 nm	4 nm	100 mm	Via PC	0.3-2.5A	0 - 9999 C	.004 @ 0.5 A	USB & RS232	Single	Yes	110 Volt	10
S2100+PE	325-1000 nm	4 nm	100 mm	Via PC	0.3-2.5A	0 - 9999 C	.004 @ 0.5 A	USB & RS232	Single	Yes	220 Volt	10



#### **Basic Visible, Series 1000**

# OFERATION INSTRUCTIONS I have not one with the 13 aris. I have a record with the 13 aris. SPECTROPHOTON SPECTROPHOTON MOSE MOSE MOSE MOSE TOTAL T

# OPERATION INSTRUCTIONS Some and was the 1st amount of the control of the control

S1000 shown with Cuvet Holder closed



View of \$1000 USB Data Port



S90-301 Optical glass Round Cuvets

#### **\$1000 Visible Spectrophotometers**

Precisely designed and ruggedly built, this easy-to-use spectrophotometer is ideal for student use in high schools, colleges and for general laboratory testing. It is suitable for general analysis and experiments such as Beer-Lambert, Absorption Spectrum, Transmittance, Chlorophyll, protein (biuret test) and more. Includes a visible wavelength spectrum chart on the instrument panel, allowing the user to make the correlation between a wavelength number to an actual color on the visible light spectrum. A versatile instrument designed to accept both round optical glass, or square 10 mm path length cuvets.

The built-in secondary filters reduce stray light and increase precision. The design is "Student-Proofed" making this a durable unit for years of use. Large digital display makes reading much easier and the new USB interface allows customers to collect and print data via a computer.

The model \$1000 now includes a USB port for data transfer to a PC. Optional Windows® based software is available for easy data collection and application expansion to Standard Curve and Kinetics. Data can easily be exported to Microsoft Excel® for further processing, analysis and documentation.

The spectrophotometer comes with 1/2" round tube holder, 10 mm square cuvet adapter, box of 12 round optical glass cuvets, Users Guide-Experiment manual and dust cover. It is pre-set for 110 Volt, but is user switchable to 220 Volt. The detachable power cord can easily be switched to a 220 Volt international plug/cord.

	\$1000
Wavelength Range	400-1000nm
Slit Width	20nm
Optical System	Single Beam, Grating System
	1200 lines/mm
Wavelength Accuracy	± 3nm
Wavelength Repeatabili	ity ± lnm
Stray Energy	< 2% @ 400nm
Photometric Range	0% - 100% T, 0 - 1.99 Abs
Photometric Accuracy	± 2%
Photometric Noise	± 1.0%
Data Port	USB
Light Source	Tungsten Halogen Lamp 6V/10W
Power Requirements	110-120V/60Hz Switchable
	220V-240V/50Hz
Instrument Dimensions	16"W x 12"D x 8"H
	406mm x 305mm x 205mm
Instrument Weight	15 lbs (7kg)

Item No.	Description
S1000	Model \$1000 Spectrophotometer, 20nm, 110V
S1000E	Model S1000E Spectrophotometer, 20nm, 220V
S1100-505	Replacement Tungsten Halogen Bulb, 6V / 10 Watt
S1000-401	Application Software for Windows
S1100-115	COD Vial Sample Holder
S90-301	Optical glass Round Cuvets, Pk/12
S90-304G	Optical glass Square Cuvets, Pk/2



# **\$1100 Series Visible Spectrophotometers**

S1100 Series offers your choice of 10 and 20nm models and is designed for use with 1/2" (13mm) round tubes or 10mm square cuvets. Optional tube holders are also available for COD, 3/4" and 1" round tubes. An adapter is supplied at no additional charge for 10mm square cuvets.

High quality silicon photodiode detector and 1200 lines/mm grating assures high performance. Auto zero function, one-touch blanking and built-in, automatic filters for easy operation. Large digital display allows quick and easy readings.


Both models in the \$1100 Series work with optional Windows® based software for easy data collection and application expansion to Standard Curve and Kinetics, Abs. and %T. Data can be easily exported to Microsoft Excel® for further processing and analysis. The software is designed for Windows based PC's, using Windows 2000, ME, or XP.

Model S1100 measures Absorbance, Transmittance and now includes a USB port; while Model S1100RS features both USB and RS-232C ports along with Absorbance, Transmittance plus Concentration (C) and Factor (F) modes.

Each Series 1100 Spectrophotometer comes complete with a manual, 12 round optical glass cuvets, a square cuvet adapter, and dust cover.

#### **CE Approved**

	S1100
Wavelength Range	335-1000nm
Slit Width	20nm
Optical System	Single Beam, Grating System
	1200 lines/mm
Wavelength Accuracy	± 2nm
Wavelength Repeatabil	lity ± lnm
Photometric Range	0% to 125% T, 0 to 2.0 Abs
Photometric Accuracy	± 2% T
Stray Light	Less than or equal to 0.5%T
	at 340 and 400nm
Data Port	USB
Light Source	Tungsten Halogen Lamp 6V/10W
Sample Compartment	Accommodates one inch round tube
	with optional holder
Power Requirements	110-120V/60Hz Switchable
	220V-240V/50Hz
<b>Instrument Dimensions</b>	16"W x 12"D x 8"H
	406mm x 305mm x 205mm
Instrument Weight	13 lbs (6 kg)







Close-up view of \$1100 instrument panel



View of S1100 USB Data Port

Se habla español. Visítenos en el Internet: www.unicosci.com/espanol o llámenos por teléfono al 609-240-5507.



#### Visible Spectro, Series 1100RS

# OFERATION INSTRUCTIONS - Foreign and was for 16 min. - Foreign and was for 16 min. - Install sample 2 - Install sample 2 - Install sample 2 - Install sample 3 - Ins



Close-up view of S1100RS instrument panel



S1100RS shown with Cuvet Holder closed



View of S1100RS USB and RS232 Data Ports

### **S1100RS Series Visible Spectrophotometers**

S1100RS Series features an accurate 10nm bandpass model and is designed for use with 1/2" (13mm) round tubes or 10mm square cuvets. Optional tube holders are also available for COD, 3/4" and 1"round tubes. An adapter is supplied at no additional charge for 10mm square optical glass cuvets.

High quality silicon photodiode detector and 1200 lines/mm grating assures high performance. Auto zero function, one-touch blanking and built-in, automatic filters for easy operation. Large digital display allows quick and easy readings.

Both models, the S1100RS and S1100RSE work with optional Windows® based software for easy data collection and application expansion to Standard Curve and Kinetics, Abs. and %T and Concentration, etc. Data can be easily exported to Microsoft Excel® for further processing and analysis. The software is designed for Windows based PC's, using Windows 2000, ME, or XP.

Models S1100RS and S1100RSE feature both USB and RS-232C ports along with Absorbance, Transmittance, Concentration (C) and Factor (F) modes. Each mode can be quickly selected with the touch of a button (please see the close-up photo of the instrument control panel on this page).

Each \$1100RS Series spectrophotometer comes complete with an Operation Manual, 12 round optical glass cuvets, a square cuvet adapter, and dust cover.

#### **CE Approved**

	S1100RS
Wavelength Range	335-1000nm
Slit Width	10nm
Optical System	Single Beam, Grating System
	1200 lines/mm
Wavelength Accuracy	± 2nm
Wavelength Repeatabil	ity ± 1nm
Photometric Range	0% to 125% T, 0 to 2.0 Abs
	0 to 1999C (0 to 1999F)
Photometric Accuracy	± 1% T
Stray Light	Less than or equal to 0.5%T
	at 340 and 400nm
Data Port	USB and RS-232
Light Source	Tungsten Halogen Lamp 6V/10W
Sample Compartment	Accommodates one inch round tube
	with optional holder
Power Requirements	220V-240V/50Hz
	110-120V/60Hz Switchable
Instrument Dimensions	16"W x 12"D x 8"H
	406mm x 305mm x 205mm
Instrument Weight	13 lbs (6 kg)

Item No.	Description
S1100RS	Model S1100RS Spectrophotometer, 10nm, 110V
S1100RSE	Model S1100RSE Spectrophotometer, 10nm, 220V
S1100-401	UNICO Application Software for PC's
S90-301	Round Optical Glass Cuvets, 10 mm, Box/12
S90-304G	Square Optical Glass Cuvets, 10 mm, Pk/2
Contact UNICO	) for additional parts and accessories



# **\$1200 Series Visible Spectrophotometers**

UNICO \$1200 Series is the best value in a precise, accurate 5nm design. High quality silicon photodiode detector and 1200 lines/mm diffraction grating assures high performance. Large digital display and built-in, automatic filters for easy operation. 1200's feature Absorbance, Transmittance, Factor and Concentration modes, also automatic zeroing and blanking with the touch of a single button.

The S1200 Series feature both USB and RS-232C ports for data transfer and you can select optional Windows® based software for easy data collection and application expansion to standard curve and kinetics, Abs. and %T. Data can be easily exported to Microsoft Excel® for further processing, analysis and storage. Software can be used with Windows 2000, ME and XP.

Bulb changes are quick and easy when needed and require no tools, or alignment. The large sample compartment can accept cuvets up to 50 mm path length and can accept a variety of optional accessories. The easy-to-use UNICO \$1200 ideal for any standard spectrophotometer application, while staying within a conservative budget.

Model \$1200 comes standard with a V-type round tube holder, a single square cuvet holder, 12 round optical glass cuvets, a set of two optical square glass cuvets, operation manual and dust cover.

Model \$1201 comes standard with a 4-cell (four position) 10 mm, square cuvet changer, set of two optical glass square cuvets, manual and dust cover.

#### **CE Approved**



\$1200 shown with Cuvet Holder closed

Item No.	Description	
S1200	Model S1200 with V-Type tube holder, 5nm, 110V	
S1200E	Model S1200E with V-Type tube holder, 5nm, 220V	
S1201	Model S1201 with 4 position cell holder, 5nm, 110V	
S1201E	Model \$1201E with 4 position cell holder, 5nm, 220V	
S1200-401	UNICO Application Software for PC's	
S90-301	Round Optical Glass Cuvets, 10 mm, Box/12	
S90-304G	Square Optical Glass Cuvets, 10 mm, Pk/2	
Contact UNICO for pricing, additional parts and accessories		




\$1200's 50mm pathlength sample chamber accepts a wide variety of cells and accessories



Close-up view of \$1200 instrument panel

	<b>S1200</b>
Wavelength Range	325-1000nm
Slit Width	5nm
Optical System	Single Beam, Grating System 1200 lines/mm
Wavelength Accuracy	± 2nm
<b>Wavelength Repeata</b>	bility ± 1nm
Photometric Range	0% T to 125% T
	0 to 2.0 A
	1 to 999 C (0 to 1999 F)
Photometric Accuracy	±1% T
Stray Light Le	ss than or equal to 0.5%T at 340 and 400nm
Light Source	Tungsten Halogen Lamp 6V/10W
Data Ports	USB and RS-232C
Sample Compartmen	Accommodates 50mm pathlength cuvet
	with optional holder
<b>Power Requirements</b>	110-120V/60Hz Switchable / 220V-240V/50H
Instrument Dimension	ns 16"W x 12"D x 7"H
	408mm x 308mm x 180mm
Instrument Weight	15 lbs (6.5 kg)





expanding on the capabilities of its analog cousin, the \$1200. The \$1205 is an all digital, fully programmable spectrophotometer with 5nm accuracy and push button, digital key pad wavelength change. This versatile instrument features a wide sample compartment and can accept cuvets up to 50 mm path-length and a variety of accessories. The wide path length and optional accessories make this unit ideal for water testing, environmental and petroleum industry testing.

The \$1205 has a high quality silicon photodiode detector and 1200 lines/mm grating assures high performance. Easy to read four (4) line digital display and built-in filters make your tests a snap. Features quick and easy automatic blanking with the push of a button.

Turn your \$1205 into a very cost effective scanning spectrophotometer with the optional Windows based software package. The \$1205 is capable of storing and recalling up to 200 standard curves and 500 test data results in its integrated memory. The \$1205's feature both a USB and R\$232 ports for interfacing with a computer. The \$1205 features optional Windows® based application software for easy data collection and application expansion ranging from wavelength scanning; to standard curve, kinetics, Abs. and %T. Data can be easily exported to Microsoft Excel® for further processing, analysis and storage. Software is designed for Windows 2000, ME and XP.

With a variety of optional accessories, the easy-touse \$1205 is ideal for any standard spectrophotometer applications, while staying within a conservative budget. Model \$1205 comes standard with a 4-cell square 10 mm path length cuvet changer, set of two optical glass square cuvets, manual and dust cover.

Wavelength Scanning Software:

Use the optional software to do a complete wavelength scan, or set the range of scan that best fits your application.

- 1. Automatically record peaks and valleys
- 2. Re-scaling axes, curves

#### **CE** approved

Item No.	Description
S1205	Model \$1205 with 4 position cell holder, 5nm, 110V
S1205E	Model \$1205E with 4 position cell holder, 5nm, 220V
S1205-401	UNICO Application Software for PC's
S90-301	Round Optical Glass Cuvets, 10 mm, Box/12
S90-304G	Square Optical Glass Cuvets, 10 mm, Pk/2





Close-up view of \$1205 instrument panel with four-line digital display.



USB and RS-232 ports.

#### **\$1205** Features and Specifications:

- Self-Calibration Feature: Instrument can perform an automated system check and calibration covering: Bulb, filter positioning, and wavelength check
- Utility Folder: From the List, the user can select Recalibrate System, Check Dark Current for diagnostic purposes, re-set time, and more
- Programmability:
  - 1. Can pre-pre-program up to 200 test methods, and store 500 results
  - 2. Program Test and Workflow Management: Can create list(s) from the "RUN TEST LIST Menu", of favorite tests, or most commonly used tests, or for a specific sequence of tests you plan to run, etc.
- USB port for data transfer and software upgrades (upgrades are no charge).
- Memorize the settings of last test, wavelength, settings, filter, etc.
- With the S-1205 Windows based software, has a scanning feature of the full wavelength range. User can define the scanning range, starting and ending wavelength, or do a full visible wavelength range scan from 325-1000 nm

	S1205
Wavelength Range	325 - 1000 nm Visible Wavelength
Slit Width	5 nm
Optical System	Single Beam, Grating System 1200 lines/mm
Wavelength Accurac	<b>y</b> ± 2nm
Wavelength Repeate	ability ± 1nm
<b>Wavelength Resoluti</b>	ion (display) Inm
Photometric Accurac	<b>y</b> ±0.004 A at 0.5 A
Photometric Range	0% T to 125% T
	-0.3 to 2.5 A
	-9999 to 9999
Photometric Repeate	<b>ability</b> ±0.003 A at 0.5 A
Stray Light	Less than or equal to 0.3% T at 340 nm
Stability	±0.002 A/hr at 500 nm after 1 hr. warm-up
Baseline flatness	±0.010 A with software
Display	LCD (4 line x 20 characters)
Light Source	Tungsten Halogen Lamp
Control Buttons	9 buttons
Printer Port	RS-232
<b>PC Communication P</b>	Port USB
Cuvet Capacity	Up to 50 mm pathlength (with longpath
	cell holder, optional)
Cell Holder	Single or 4 cell, manual
Abs / %T Mode	Display: Abs and % T
	Standard curve: 1st through zero, 1st, Factor
Memory Capability	Up to 200 standard curves
Utility	Dark current, system self calibration
<b>Power Requirements</b>	
<b>Instrument Dimensio</b>	ons 20"W x 17"D x 8"H
	420mm x 340mm x 180mm
Instrument Weight	19 lbs (8.5 kg)
Warranty	One year

### Accessories for Models \$1200 and \$1205



#### Universal Test Tube Holder (Cat. No. 51200-101)

Universal test tube holder with base for 8 to 25mm diameter test tubes including COD.

S1200-101 with cell

#### Single Cell Holder (Cat. No. \$1200-102)

Single cell holder with base for 10mm square cuvets.



\$1200-102



#### Four-Position Cell Holder (Cat. No. S1200-103)

Four-position cell holder for 10mm square cuvets.



#### Longpath Cell Holder (Cat. No. 51200-104)

Longpath four-position cell holder for up to 50mm rectangular cuvets.



\$1200-104



#### Water-jacketed Cell Holder (Cat. No. \$1200-105)

Water-jacketed 10mm single cell holder with base and panel for temperature control applications (requires water bath).

S1200-105

#### Two-in-One Cell Holder (Cat. No. \$1200-108)

Holds either a 100mm square or a 10mm round cuvet.



S1200-108









Close-up view of optional built-in printer.

	S2100+ / S2100+P
Wavelength Range	325-1000 nm Visible Wavelength
Slit Width	4 nm
Optical System S	ingle Beam, Grating System 1200 lines/mm
Wavelength Accuracy	± 2nm
Wavelength Repeatal	bility ± 1nm
Wavelength Resolution	on (display) Inm
Photometric Accuracy	±0.004 A at 0.5 A
Photometric Range	0% T to 125% T
	-0.3 to 2.5 A
	-9999 to 9999
Photometric Repeatal	<b>bility</b> ±0.003 A at 0.5 A
Stray Light	Less than or equal to 0.3% T at 340 nm
Stability ±	-0.002 A/hr at 500 nm after 1 hr. warm-up
Baseline flatness	±0.010 A with PC software
Display	LCD (4 line x 20 characters)
Control Buttons	9 buttons, Key Pad
Printer (optional)	Built-in Thermal Printer (optional)
Communication Ports	USB and RS-232
Cuvet Capacity	Up to 100 mm pathlength (with longpath
	cell holder, optional)
Cell Holder	Single or 4 cell, manual
Abs / %T Mode	Display: Abs and % T
	tandard curve: 1st through zero, 1st, Factor
Memory Capability	Up to 200 standard curves
	ark current, D2 on/off, system self calibration
Power Requirements	100-240 VAC, 100W
Instrument Dimension	
	420mm x 340mm x 180mm
Instrument Weight	22 lbs (10 kg)
Warranty	One year

#### S2100+ and S2100+P Visible Programmable Spectrophotometers

UNICO is proud to introduce the new S2100+ Series spectrophotometers. This new series offers software-based wavelength scanning and models with built-in printers (S2100+P). These new spectros are programmable, affordable, and simple to use; while providing excellent results.

The large 4-line LCD display can be read from any angle. The touch buttons make commanding easy and data entry quick and convenient. The menu is task driven and makes using your 2100+ user-friendly. A large on-board memory is capable of storing up to 200 test methods (standard curves) and saving up to 500 test data.

Each unit comes standard with both USB and RS-232 data ports for communication. Two optional features make documentation and reporting a snap; optional software for PC's which includes scanning capability, and units with Built-in printers (S2100+P series).

The USB port can be used for both data transfer and software upgrades. Periodically we make system upgrades, and you will be able to update the functionality of your S2100+ at no charge via the USB port whenever upgrades are available; please contact UNICO for details.

The optional Windows® XP based software expands the capabilities of the instrument including Standard Curves, Kinetics, Abs. & %T, and wavelength scanning. Data can be easily exported to Excel® for further processing, analysis and storage.

The S2100+ comes standard with an operation manual, dust cover, and a set of 4 10 mm square optical glass cuvets.

#### **CE Approved**



Close-up view of \$2100+ instrument panel with four-line digital display.

Item No.	Description
S2100+	Model S2100+ with 4 position cell holder, 4nm, 110V
S2100+P	Model S2100+P with Built-in Printer, 4nm, 110V
S2100+E	Model S2100+E with 4 position cell holder, 4nm, 220V
S2100+PE	Model S2100+PE with Built-in Printer, 4nm, 220V
S2100-401	UNICO Application Software with WL Scanning for PC's
S90-304G	Square Optical Glass Cuvets, 10 mm, Pk/2
S90-309Q	Square Quartz UV Transparent Cuvets, 10 mm, Pk/2
Contact UNICC	of for additional parts and accessories

See pages 13 and 14 for S2100+ software and accessories.



#### **UV-Visible Spectro Selection Guide**

#### **UV-Visible Spectrophotometer Selection Guide**



S2100UV+

#### **Electromagnetic Spectrum:**

Wavelength Range in Nanometers (nm):

Ultra-Violet (UV), Visible (Vis) and Infrared Light (IR)

Ultra-Violet (UV) Spectrum: 200nm – 400nm
 Visible Light Spectrum (Vis): 400nm – 750nm

• Infrared Spectrum (IR): 750nm - 300,000nm



SQ3802



SQ2800

#### **UNICO Visible Light Spectrophotometers**

Wavelength Range by series:

\$2100UV+ Series: 200nm - 1,000nm
 \$Q2800 Series: 190nm - 1,100nm
 \$Q3800 Series: 190nm - 1,100nm
 \$Q4800 Series: 190nm - 1,100nm



SQ4802

How do I select the right spectrophotometer for my applications? We can help. Contact UNICO by phone, fax, or e-mail and our technicians will be glad to help select the right instrument for your needs. Please see the Selection Guide below. This will help you understand the general specifications by model.

Model	Wavelength Range (nm)	Slit Width Bandpass	Cell Pathlength		Absorbance Range	Concentration Conc., C	Photometric Accuracy	Data Port(s)	Light Beam	Printer	Voltage Pre-set	Pg.
S2100UV+	200-1000 nm	4 nm	100 mm	Via PC	0.3-2.5A	0 - 9999	.004 @ 0.5 A	USB & RS232	Single	Optional	110 Volt	12
S2100UV+E	200-1000 nm	4 nm	100 mm	Via PC	0.3-2.5A	0 - 9999	.004 @ 0.5 A	USB & RS232	Single	Optional	220 Volt	12
S2100UV+P	200-1000 nm	4 nm	100 mm	Via PC	0.3-2.5A	0 - 9999	.004 @ 0.5 A	USB & RS232	Single	Built-in	110 Volt	13
52100UV+PE		4 nm	100 mm	Via PC	0.3-2.5A	0 - 9999	.004 @ 0.5 A	USB & RS232	Single	Built-in	220 Volt	13
3210001111	200 1000 11111	7 11111	100 11111	11410	0.0 2.071	0 ,,,,	.004 & 0.071	OOD O NOZOZ	omgic	Dom m	220 1011	
SpectroQuest	тм											
SQ2800	190-1100 nm	4 nm	100 mm	Yes	0.3-3A	0 - 9999	.004 @ 0.5 A	USB & RS232	Single	Optional	110 Volt	16
SQ2800E	190-1100 nm	4 nm	100 mm	Yes	0.3-3A	0 - 9999	.004 @ 0.5 A	USB & RS232	Single	Optional	220 Volt	16
SQ2800P	190-1100 nm	4 nm	100 mm	Yes	0.3-3A	0 - 9999	.004 @ 0.5 A	USB & RS232	Single	Built-in	110 Volt	16
SQ2800PE	190-1100 nm	4 nm	100 mm	Yes	0.3-3A	0 - 9999	.004 @ 0.5 A	USB & RS232	Single	Built-in	220 Volt	16
SQ2802	190-1100 nm	1.8 nm	100 mm	Yes	0.3-3A	0 - 9999	.004 @ 0.5 A	USB & RS232	Single	Optional	110 Volt	1 <i>7</i>
SQ2802E	190-1100 nm	1.8 nm	100 mm	Yes	0.3-3A	0 - 9999	.004 @ 0.5 A	USB & RS232	Single	Optional	220 Volt	17
SQ2802S	190-1100 nm			Yes	0.3-3A	0 - 9999	.004 @ 0.5 A	USB & RS232	Single	Optional	110 Volt	17
SQ2802SE	190-1100 nm	1, 1.8, 4 nm	100 mm	Yes	0.3-3A	0 - 9999	.004 @ 0.5 A	USB & RS232	Single	Optional	220 Volt	17
SQ2802PCS	190-1100 nm	1 1 0 4	100 mm	Yes	0.3-3A	0 - 9999	.004 @ 0.5 A	USB & RS232	C:nala	Optional	110 Volt	17
	190-1100 nm	, ,	100 mm	Yes	0.3-3A 0.3-3A	0 - 9999	.004 @ 0.5 A	USB & RS232	Single	Optional	220 Volt	
3Q26U2PC3E	190-1100 nm	1, 1.0, 4 nm	100 mm	res	0.3-3A	0 - 9999	.004 @ 0.5 A	USB & RSZSZ	Single	Opilonal	220 VOII	17
SQ3802	190-1100 nm	1.8 nm	100 mm	Yes	0.3-3A	0 - 9999	.004 @ 0.5 A	USB & RS232	Split	Optional	110 Volt	18
SQ3802E	190-1100 nm	1.8 nm	100 mm	Yes	0.3-3A	0 - 9999	.004 @ 0.5 A	USB & RS232	Split	Optional	220 Volt	18
SQ4802	190-1100 nm	1.8 nm	100 mm	Yes	0.3-3A	0 - 9999	.004 @ 0.5 A	USB & RS232	Double	Optional	110 Volt	19
SQ4802E	190-1100 nm	1.8 nm	100 mm	Yes	0.3-3A	0 - 9999	.004 @ 0.5 A	USB & RS232	Double	Optional	220 Volt	19









Close-up view of optional built-in printer.

	S2100UV+
Wavelength Range	200 - 1000 nm Visible Wavelength
Slit Width	4 nm
Optical System Si	ingle Beam, Grating System 1200 lines/mm
Wavelength Accuracy	± 2nm
Wavelength Repeatak	bility ± 1 nm
Wavelength Resolutio	n (display) 1 nm
<b>Photometric Accuracy</b>	±0.004 A at 0.5 A
Photometric Range	0% T to 125% T
	-0.3 to 2.5 A
	-9999 to 9999
Photometric Repeatab	bility ±0.003 A at 0.5 A
Stray Light	Less than or equal to 0.3% T at
	220 nm and 340 nm
Stability ±	0.002 A/hr at 500 nm after 1 hr. warm-up
Baseline flatness	±0.010 A with PC software
Display	LCD (4 line x 20 characters)
Control Buttons	9 buttons, Key Pad
Printer (optional)	Built-in Thermal Printer (optional)
<b>Communication Ports</b>	USB and RS-232
Cuvet Capacity	Up to 100 mm pathlength (with longpath
	cell holder, optional)
Cell Holder	Single or 4 cell, manual
Abs / %T Mode	Display: Abs and % T
<b>Quantitative Mode</b> S	tandard curve: 1st through zero, 1st, Factor
Memory Capability	Up to 200 standard curves
<b>Utility</b> Do	ark current, D2 on/off, system self calibration
Power Requirements	100/220V switchable, 150W
<b>Instrument Dimension</b>	20"W x 17"D x 8"H
	420mm x 340mm x 180mm
Instrument Weight	22 lbs (10 kg)
Warranty	One year

# **S2100UV+ UV/Vis Spectrophotometers**

UNICO is proud to introduce the new \$2100UV+ Series UV/Vis spectrophotometers. This new series offers software-based wavelength scanning and models with built-in printers (\$2100UV+P). These new spectros are programmable, affordable, and simple to use; while providing excellent results.

The large 4-line LCD display can be read from any angle. The touch buttons make commanding easy and data entry quick and convenient. The menu is task driven and makes using your 2100UV+ user-friendly. A large on-board memory is capable of storing up to 200 test methods (standard curves) and saving up to 500 test data.

Each unit comes standard with both USB and RS-232 data ports for communication. Two optional features make documentation and reporting a snap; optional software for PC's which includes scanning capability, and units with Built-in printers (\$2100+P series).

The USB port can be used for both data transfer and software upgrades. Periodically we make system upgrades, and you will be able to update the functionality of your S2100+ at no charge via the USB port whenever upgrades are available; please contact UNICO for details.

The optional Windows® XP based software expands the capabilities of the instrument including Standard Curves, Kinetics, Abs. & %T, DNA/Protein ratio and wavelength scanning. Data can be easily exported to Excel® for further processing, analysis and storage.

The S2100UV+ comes standard with an operation manual, dust cover, a set of two 10 mm UV Transparent Quartz cuvets and a set of four 10 mm square optical glass cuvets.

#### **CE Approved**



Close-up view of S2100UV+ instrument panel with four-line digital display.

Item No.	Description
S2100UV+	Model S2100UV+ with 4 position cell holder, 4nm, 110V
S2100UV+P	Model S2100UV+P with Built-in Printer, 4nm, 110V
S2100UV+E	Model S2100UV+E with 4 cell holder, 4nm, 220V
S2100UV+PE	Model S2100UV+PE with Built-in Printer, 4nm, 220V
S2100-401	UNICO Application Software with WL Scanning for PC's
S90-304G	Square Optical Glass Cuvets, 10 mm, Pk/2
S90-309Q	Square Quartz UV Transparent Cuvets, 10 mm, Pk/2
Contact UNICO	for additional parts and accessories



# **S2100+ Series**Options and Accessories

#### Test Tube Holder (Item #S2100-101P)

Test tube holder kit for 8-20mm diameter test tubes. Includes universal base, V-type tube holder. The maximum tube height is 100mm.



S2100-101P

S2100-102P

#### Long Path Cell Holder (Item #S2100-102P)

Rectangular long path cell holder kit for single cell up to 100mm pathlength. Includes a universal base and one cell holder.

#### Single Square Cell Holder (Item #S2100P-103P)

Single square cell holder 10mm pathlength cuvet.



S2100P-103P

S2100-104P with cell

#### Cylindrical Cell Holder (Item #S2100-104P)

Cylindrical cell holder kit for single cell up to 100mm pathlength (20mm diam.). Includes a universal base and one holder.

#### Water-Jacketed Cell Holder (Item #S2100P-105P)

Water-Jacketed single cell holder kit including a universal base and one water-jacketed cell holder for 10mm square cuvet. It maintains desired temperature by circulating constant-temperature water from a Water Bath (water bath required but not included).



S2100P-105P



S2100-106P

#### Micro Cell Holder (Item #S2100-106P)

Measure a sample with small volumes of 100 uL using micro cell holder. The x-y adjustable mechanism is used to align cell with optical beam for optimized results.

#### Peltier Unit (Item #S2100-107P)

Peltier unit for continuous temperature control from 15° to 40° C. The x-y adjustable mechanism is used to align cell with optical beam for micro flow cell setup. The temperature display resolution is 0.1° C. The unit consists of a controller and a thermoelectrically controlled cell holder.



S2100-107P

S2100P-108P

#### Ambient Sipper Unit (Item #S2100P-108P)

Sipper system for single cell flow thru. The x-y adjustable mechanism is used to align cell with optical beam for micro flow cell setup. The sipper unit consists of a flow-thru controller with peristaltic pump and a flow-thru front panel (flow cell and tubing not included). Note: requires flow cell and proper tubing to complete flow thru setup.

#### Peltier/Sipper System (Item#S2100P-109P)

Peltier/Sipper system for single 10mm cell flow thru and continuous temperature control from 15° to 40° C. The x-y adjustable mechanism is used to align cell with optical beam for micro flow cell setup. The



S2100P-109P

temperature display resolution is 0.1° C. The unit consists of a Petier/Sipper controller with peristaltic pump and a thermoelectrically temperature controlled cell holder with panel. The unit can be used as flow thru only or temperature control only. Note: requires flow cell and proper tubing to complete flow thru setup.



S2100-110P

#### Reflectance Measurement Attachment (5° incident angle) (Item #S2100-110)

The technique of reflectance measurement is used for evaluation of materials relative to a reflectance surface. The minimum sample is (L) 30 x (W) 30 mm.

#### Four-Position Cell Holder (Item #S2100-111)

Four-position cell holder for 10mm square cuvets.



S2100-111



#### S2100+ and S2100UV+ Series Spectrophotometer Application Software (S2100P-401)


This software is designed for Windows® 2000, ME, or XP operating systems. Install the Windows based software and start collecting your data.

The UNICO Application Software expands your applications, assists in your data documentation and provides complete control of the spectrophotometer from a computer, rather than the instrument panel.

Applications include:


- 1. Abs./%T/Conc.
- 2. Standard Curve (Quantitative)
- 3. Kinetics (Abs. vs. Time)
- 4. Scanning
- 5. DNA/Protein Ratio (UV Range only)

#### Abs/%T/Conc



Measure absorbance, transmittance or concentration with standard or known factor.


#### **Quantitative (Standard Curve)**



Use up to 8 standards to establish standard curves. Features four methods for fitting a curve:

- 1. Linear fit
- 2. Linear fit through zero
- 3. Square fit
- 4. Segmented


#### **Wavelength Scanning Software**




- 1. Full or partial wavelength scans; set the range required by application
- 2. Automatically record peaks and valleys
- 3. Re-scaling axes, curves

#### **Kinetics (Abs vs. Time)**

The Kinetics mode may be used for time course scanning or reaction rate calculations. Abs. vs. time graphs is displayed on the screen in real time and sayed.



Wait time, measurement time and time intervals may be entered.



#### **DNA/Protein**

Concentration and DNA purity are calculated: Absorbance ratios 260nm/280nm with optional subtracted absorbance at 320nm.

DNA Concentration = 62.9 x A260 - 36.0 x A280

Protein Concentration = 1552 x A260 - 757.3 x A280

Additional wavelengths and factors may be programmed and saved. Please note: this feature can only be used with S2100UV+ Series spectrophotometers.



#### **SpectroQuest**<sup>TM</sup> **Line**

#### SpectroQuest™ Line of UV/Vis Scanning Spectrophotometers

UNICO is proud to introduce the all-new complete SpectroQuest™ line of UV-Visible spectrophotometers. The SpectroQuest line consists of four distinct series to meet the broad requirements of education, industrial and research applications:

- SQ-2800 Series Single Beam 4nm UV/Vis spectrophotometer
- SQ-2802 Series Single Beam 1.8nm or variable slits UV/Vis spectrophotometer
- SQ-3802 Series Split Beam 1.8nm UV/Vis spectrophotometer
- SQ-4802 Series Double Beam 1.8nm UV/Vis spectrophotometer

All SpectroQuest spectrophotometers feature high performance sealed optics mounted on a stable machined platform. The innovative optical layout and state of art monochromator with high-grade blazed holographic grating ensure accuracy. Its integrated design assures long-term stability and durability. The precisely aligned detector and quality deuterium and halogen lamps enhance the precision across the UV/Vis spectrum starting from 190nm and into the near-infrared (NIR) 1100nm. The comprehensive features, sophisticated powerful software, variety of accessories and model configurations will meet or exceed your expectations for performance and value.



- Choice of single beam, split beam or double beam designs
- Fixed or variable slits (bandwidths)
- PC models or Stand-alone models with large LCD display
- Optional built-in printer on SQ2800 models
- Non-volatile memory storage and one-button easy recall
- Sealed keypad with alpha-numeric entry for user file names and settings
- Pre-aligned deuterium lamp for easy lamp replacement. Lamp usage and the status of the lamps may be monitored
- All models with USB and RS-232 ports for communication and software upgrades
- Powerful built-in or PC Windows software including sophisticated utility programs
- Data Download-to-PC software for stand-alone models (optional)
- Software upgrade capability
- Real-time clock for date and time stamping of results
- Performance validation and report (GLP compliance)
- Full CE compliance

UNICO continues to use Tungsten-Halogen (Visible) and Deuterium (UV) light sources for the added stability these light sources provide as compared to Xenon sources.

SpectroQuest<sup>™</sup> is a registered trademark of United Products and Instruments, Inc., dba UNICO.



SQ-2800



SQ-2802



SQ-3802



SQ-4802



#### SpectroQuest™ Series 2800





SQ-2800 shown with 100 mm sample compartment open

	SQ-2800
Wavelength Range	190 - 1100 nm
Slit Width	4 nm
Optical System Si	ngle Beam, Grating System 1200 lines/mm
Wavelength Accuracy	± 0.8 nm
<b>Wavelength Resolution</b>	± 0.1 nm
Wavelength Repeatabil	lity ± 0.5 nm
Photometric Range	0% T to 200% T
	-0.3 to 3 A
	0 to 9999 Conc
Photometric Accuracy	± 0.5% T
Photometric Repeatabil	lity Better than 0.3% T
Stray Light	Less than 0.15% T
Baseline Flatness	± 0.004 A
Stability	0.002 A/h @ 500 nm
Scanning Speed	High, Medium, Low
	Maximum 1000 nm/minute
Light Source	Halogen, Deuterium (pre-aligned)
Display	Graphic LCD (320 x 240) dots
Keyboard Control	29 Membrane keypad
Data Output	USB and RS232C, Parallel printer port
Sample Compartment	Accommodates 100mm pathlength
	cuvet with optional holder
Printer	Mini parallel printer optional
Power Requirements	110-120V/60Hz Switchable
	220V-240V/50Hz
Instrument Dimensions	21.7"W x 16.5"D x 10.6"H
	550mm x 420mm x 270mm
Instrument Weight	44 lbs (20 kg)
Warranty	One Year

# SQ-2800 Single Beam UV/Vis Spectrophotometers

SQ-2800 is the most economic general-purpose design in the SpectroQuest line. It is a stand-alone model with 4nm fixed bandwidth and has all the features that SpectroQuest line offers for a stand-alone unit. It provides excellent performance for measurements in the range of 190nm to 1100nm. It has a large angled LCD screen with contrast adjustment for comfortable viewing in a variety of room conditions. The large sample compartment (100 mm) accommodates a wide range of cell holders and accessories including peristaltic sipper and peltier system. Optional PC download software and PC Windows® application software make this instrument very versatile. SQ2800P Models are available with built-in printers.

	SQ-2800P
Wavelength Range	190 - 1100 nm
Slit Width	4 nm
Optical System Si	ingle Beam, Grating System 1200 lines/mm
Wavelength Accuracy	± 0.8 nm
<b>Wavelength Resolution</b>	
Wavelength Repeatabi	<b>lity</b> ± 0.5 nm
Photometric Range	0% T to 200% T
-	-0.3 to 3 A
	0 to 9999 Conc
Photometric Accuracy	± 0.5% T
Photometric Repeatabil	lity Better than 0.3% T
Stray Light	Less than 0.15% T
Baseline Flatness	± 0.004 A
Stability	0.002 A/h @ 500 nm
Scanning Speed	High, Medium, Low
	Maximum 1000 nm/minute
Light Source	Halogen, Deuterium (pre-aligned)
Display	Graphic LCD (320 x 240) dots
Keyboard Control	29 Membrane keypad
Data Output	USB and RS232C, Parallel printer port
Sample Compartment	Accommodates 100mm pathlength
	cuvet with optional holder
Printer	Built-in
Power Requirements	110-120V/60Hz Switchable
	220V-240V/50Hz
<b>Instrument Dimensions</b>	21.7"W x 16.5"D x 10.6"H
	550mm x 420mm x 270mm
Instrument Weight	44 lbs (20 kg)
Warranty	One Year

Item No.	Description
SQ2800	Model SQ2800 Scanning UV/Vis, 4nm, 110V
SQ2800E	Model SQ2800E Scanning UV/Vis, 4nm, 220V
SQ2800P	Model SQ2800P with Built-in Printer, 4nm, 110V
S2Q2800PE	Model SQ2800PE with Built-in Printer, 4nm, 220V
SQ2800-405	SpectroQuest Data Download Software for PC's
SQ2800-401	Advanced Application Software for PC's
Contact UNICO	for additional parts and accessories

Did you know that UNICO services spectrophotometers? Contact us for more details.



# **SQ-2802 Series Single Beam UV/Vis Spectrophotometers**

SQ-2802 series is an advanced single beam design consisting of three models: Stand-alone model SQ-2802 with 1.8nm fixed bandpass and model SQ-2802S with variable slits (1nm, 1.8nm, and 4nm); PC model SQ-2802PC with 1.8nm fixed bandpass and can be controlled either by the keypad, or remote controlled from a computer with the included software package.

SQ-2802 has all the features that SpectroQuest line offers in a stand-alone unit. The PC models come standard with Windows® based application software (PC not included) and can perform basic Abs./%T/Conc. tests without a PC. All instruments provide excellent performance and flexibility for your applications. The large sample compartment (100 mm) accommodates a wide range of cell holders and accessories including peristaltic sipper and peltier system.



	SQ-2802 / SQ-2802S
Wavelength Range	190 - 1100 nm
Slit Width	1.8 nm (SQ-2802)
	1, 1.8, 4 nm (SQ-2802S)
Optical System Sir	ngle Beam, Grating System 1200 lines/mm
Wavelength Accuracy	± 0.3 nm
<b>Wavelength Resolution</b>	
<b>Wavelength Repeatab</b>	ility ± 0.2 nm
Photometric Range	0% T to 200% T
	-0.3 to 3 A
	0 to 9999 Conc
Photometric Accuracy	± 0.3% T
Photometric Repeatabi	<b>ility</b> Better than 0.2% T
Stray Light	Less than 0.10% T
Baseline Flatness	± 0.002 A
Stability	0.002 A/h @ 500 nm
Scanning Speed	High, Medium, Low
	Maximum 1000 nm/minute
Light Source	Halogen, Deuterium (pre-aligned)
Display	Graphic LCD (320 x 240) dots
Keyboard Control	29 Membrane keypad
Data Output	USB and RS232C, Parallel printer port
Sample Compartment	Accommodates 100mm pathlength
	cuvet with optional holder
Printer	Mini parallel printer optional
Power Requirements	110-120V/60Hz Switchable
	220V-240V/50Hz
Instrument Dimensions	24.5"W x 15.7"D x 11"H
	620mm x 400mm x 280mm
Instrument Weight	48 lbs (22 kg)
Warranty	One Year

Item No.	Description
SQ2802	Model SQ2802 Scanning UV/Vis, 1.8nm, 110V
SQ2802E	Model SQ2802E Scanning UV/Vis, 1.8nm, 220V
SQ2802S	Model SQ2802S Scanning UV/Vis, Variable nm, 110V
S2Q2802SE	Model SQ2802SE Scanning UV/Vis, Variable nm, 220V
SQ2800-405	SpectroQuest Data Download Software for PC's
SQ2800-401	Advanced Application Software for PC's
Contact UNICO	for additional parts and accessories

	SQ-2802PC / SQ-2802PCS
Wavelength Range	190 - 1100 nm
Slit Width	1.8 nm (SQ-2802PC)
	1, 1.8, 4 nm (SQ-2802PCS)
	ingle Beam, Grating System 1200 lines/mm
<b>Wavelength Accuracy</b>	± 0.3 nm
Wavelength Resolution	
<b>Wavelength Repeata</b>	
Photometric Range	0% T to 200% T
	-0.3 to 3 A
	0 to 9999 Conc
Photometric Accuracy	
Photometric Repeatal	bility Better than 0.2% T
Stray Light	Less than 0.10% T
Baseline Flatness	± 0.002 A
Stability	0.002 A/h @ 500 nm
Scanning Speed	600 nm/minute
Light Source	Halogen, Deuterium (pre-aligned)
Display	20 x 2 LCD
Keyboard Control	8 Key Buttons (inactive when PC connected)
Data Output	USB and RS232C, Parallel printer port
Sample Compartmen	
	cuvet with optional holder
Printer	Mini parallel printer optional
<b>Power Requirements</b>	220V-240V/50Hz
	110-120V/60Hz Switchable
Instrument Dimension	15 24.5"W x 15.7"D x 11"H
	620mm x 400mm x 280mm
Instrument Weight	53 lbs (24 kg)
Warranty	One Year

Item No.	Description
SQ2802PC	Model SQ2802PC Scanning UV/Vis, 1.8nm, 110V
SQ2802PCE	Model SQ2802PCE Scanning UV/Vis, 1.8nm, 220V
SQ2802PCS	Model SQ2802PCS Scanning UV/Vis, Variable nm, 110V
SQ2802PCSE	Model SQ2802PCSE Scanning UV/Vis, Variable nm, 220V
SQ2800-405	SpectroQuest Data Download Software for PC's
SQ2800-401	Advanced Application Software for PC's
Contact UNICO	for additional parts and accessories



# SQ-3802 Split Beam UV/Vis Spectrophotometers



SQ-3802 is a split beam, scanning design. It is a stand-alone model with 1.8nm fixed bandwidth and has all the features that SpectroQuest line offers in a stand-alone unit. The second detector is simultaneously monitoring the system stability to optimize measurement accuracy. It provides excellent performance for measurements in the range of 190nm to 1100nm. It has a large angled LCD screen with contrast adjustment. The large sample compartment (100 mm) accommodates a wide range of cell holders and accessories including peristaltic sipper and peltier system. Optional PC download software and PC Windows® application software make this a versatile instrument.

	SQ-3802
Wavelength Range	190 - 1100 nm
Slit Width	1.8 nm
Optical System	Split Beam, Grating System 1200 lines/mm
<b>Wavelength Accuracy</b>	± 0.3 nm
Wavelength Resolution	<b>en</b> ± 0.1 nm
<b>Wavelength Repeatal</b>	
Photometric Range	0% T to 200% T
	-0.3 to 3 A
	0 to 9999 Conc
Photometric Accuracy	
Photometric Repeatal	
Stray Light	Less than 0.10% T
Baseline Flatness	± 0.002 A
Stability	0.001 A/h @ 500 nm
Scanning Speed	High, Medium, Low
	Maximum 1000 nm/minute
Light Source	Halogen, Deuterium (pre-aligned)
Display	Graphic LCD (320 x 240) dots
Keyboard Control	29 Membrane keypad
Data Output	USB and RS232C, Parallel printer port
Sample Compartment	
	cuvet with optional holder
Printer	Mini parallel printer optional
Power Requirements	110-120V/60Hz Switchable
	220V-240V/50Hz
Instrument Dimension	
	620mm x 400mm x 280mm
Instrument Weight	50 lbs (23 kg)
Warranty	One Year

Item No.	Description
SQ3802	Model SQ3802 Split Beam UV/Vis, 1.8nm, 110V
SQ3802E	Model SQ3802E Split Beam UV/Vis, 1.8nm, 220V
SQ3802-401	Advanced Application Software for PC's
Contact UNICO	for additional parts and accessories

#### Peltier/Sipper System (Cat. No. SQ2800-109P)

Peltier/Sipper system for a single cell flow thru and continuous temperature control from 15° to 40° C. The x-y adjustable mechanism is used to align cell with optical beam for micro flow cell setup. The temperature display resolution is 0.1° C. The unit consists of a Peltier/Sipper controlled with peristaltic pump and a thermoelectrically temperature controlled cell holder with panel. The unit can be used a flow thru only or temperature control only.

Note: requires flow cell and proper tubing to complete flow thru setup.



SQ2800-109P

When using a spectrophotometer in the UV range, you must use UV-transparent cuvets such as UNICO Quartz cuvets, or UV transparent plastic cuvets.

See page 26 for an extensive selection of UNICO cuvets.



SQ-4802 Double Beam UV/Vis Spectrophotometers

SQ-4802 is a double beam, scanning design. It is a stand-alone model with 1.8nm fixed bandwidth and has all the features that SpectroQuest line offers in a stand-alone unit. The two detectors are measuring both the test sample cell and reference sample cell simultaneously for optimizing measurement accuracy and stability. It provides excellent performance for measurements in the range of 190nm to 1100nm. It is suitable for pharmaceutical, biochemical and clinical lab applications as well as routine applications such as quantitative analyses, kinetics, spectrum scanning, multiple components and DNA/Protein. Optional PC download software and PC Windows® application software make this instrument versatile.



	SQ-4802			
Wavelength Range	190 - 1100 nm			
Slit Width	1.8 nm			
Optical System Do	uble Beam, Grating System 1200 lines/mm			
Wavelength Accuracy	± 0.3 nm			
<b>Wavelength Resolution</b>	± 0.1 nm			
Wavelength Repeatable	ility ± 0.2 nm			
Photometric Range	0% T to 200% T			
	-0.3 to 3 A			
	0 to 9999 Conc			
Photometric Accuracy	± 0.3% T			
Photometric Repeatabi	ility Better than 0.2% T			
Stray Light	Less than 0.10% T			
Baseline Flatness	± 0.002 A			
Stability	0.001 A/h @ 500 nm			
Scanning Speed	High, Medium, Low			
	Maximum 1000 nm/minute			
Light Source	Halogen, Deuterium (pre-aligned)			
Display	Graphic LCD (320 x 240) dots			
Keyboard Control	29 Membrane keypad			
Data Output	USB and RS232C, Parallel printer port			
Sample Compartment	Accommodates 100mm pathlength			
	cuvet with optional holder			
Printer	Mini parallel printer optional			
Power Requirements	110-120V/60Hz Switchable			
	nts 110-120V/60Hz Switchable 220V-240V/50Hz			
Instrument Dimensions	24.5"W x 15.7"D x 11"H			
	620mm x 400mm x 280mm			
Instrument Weight	± 0.1 nm  by ± 0.2 nm  0% T to 200% T  -0.3 to 3 A  0 to 9999 Conc  ± 0.3% T  by Better than 0.2% T  Less than 0.10% T  ± 0.002 A  0.001 A/h @ 500 nm  High, Medium, Low  Maximum 1000 nm/minute  Halogen, Deuterium (pre-aligned)  Graphic LCD (320 x 240) dots  29 Membrane keypad  ISB and RS232C, Parallel printer port  Accommodates 100mm pathlength  cuvet with optional holder  Mini parallel printer optional  110-120V/60Hz Switchable  220V-240V/50Hz  24.5"W x 15.7"D x 11"H			
Warranty	One Year			

Item No.	Description			
SQ4802	Model SQ4802 Double Beam UV/Vis, 1.8nm, 110V			
SQ4802E	Model SQ4802E Double Beam UV/Vis, 1.8nm, 220V			
SQ4802-401	Advanced Application Software for PC's			
	Six-Position Auto Cell Changer			
Contact UNICO for additional parts and accessories				



SQ-4802 shown with sample compartment open



SQ4802-120 Six-Position Auto Cell Holder



#### Powerful Integrated Software for Data Acquisition



#### 1. Basic Mode

Absorbance, %T Transmittance or Concentration measurements.

#### 2. Quantitative

Establish or use stored calibration equation to measure the concentration of unknowns.

#### 3. WL scan

Spectrum scan of sample at any selected wavelength range with choice of scanning speed and wavelength interval. You can also select the scan speed: Low, Medium, or High.

#### 4. Kinetics

Measurement of absorbance changing vs. time with reaction rate calculation function.

#### 5. DNA/Protein

Calculation of concentration and DNA purity. Ratio at other wavelengths can be measured.

#### 6. Multi Wavelength

Measurement at multiple wavelengths to analyze and determine the composition of the mixtures.

#### 7. Utility - GLP

Utility programs offer wavelength and photometric accuracy validation for GLP compliance. It contains useful programs and tools such as re-set dark current, re-set lamp change-over wavelength, lamp usage set, set clock, etc.

#### 8. Defined Test

This is an open platform for use defined programs, multiple test protocols can be stored.


#### **Basic Photometric Mode**

Measures Absorbance, %T and Concentration with entry of Concentration Factor or the Concentration of the standard. Units such as ug/mL, mg/mL, g/L, ppb, ppm, %, I.U., mM/L, M/L may be selected or other units may be entered via the keypad. Continuous display of the result means there is no need to press a button to read.




	Basic Mode Test Report
Wavelength	650.0nm.
Result:	0.000Abs
Date and Time:	12.05.2008 12:41:12





#### **Quantitative**

Up to 10 standard solutions may be used to establish calibration equation curve. Choice of four methods for fitting a curve through the calibration points: Linear fit, Linear fit through zero, square fit and cubic fit.





There are three kinds of correction methods:


- 1. Single wavelength method
- 2. Iso-absorbance (two wavelength method):
  The absorbance at the measurement (peak) wavelength is
  measured relative to the absorbance at a second (valley)
  wavelength. This minimizes the effects of cell difference
  and turbidity
- 3. Three-point:

  The absorbance of the peak itself is measured by subtracting the calculated tangent joining the valleys on each side of the peak



#### **SpectroQuest**<sup>TM</sup> **Operating Software**





#### **Wavelength Scanning**

The wavelength scan intervals are 0.1, 0.2, 0.5, 1, 2, 5nm, and Hi, Medium and Low scan speeds are available. Scan speeds vary from 100 to 1000 nm/min. Wavelengths are scanned from high to low so that the instrument waits at high wavelengths. This minimizes the degradation of UV sensitive samples. Precise control of filter and lamp changes means that their effects are not seen on the final scan. Post-run manipulation includes re-scaling axes, curve tracking and peak picking.



#### **DNA/Protein**

Concentration and DNA purity are calculated:


Absorbance ratios 260nm/280nm or 260nm/230nm


With optional subtracted absorbance at 320nm

DNA Concentration = 62.9 x A260 - 36.0 x A280 or 49.1 x A260 - 3.48 x A230 Protein Concentration = 1552 x A260 - 757.3 x A280 or 183 x A260 - 75.8 x A230

#### Kinetics


This mode may be used for time course scanning or reaction rate calculations. Abs. vs. time graphs are displayed on the screen in real time. Wait time and measurement time up to 12 hours may be entered with time intervals of 0.5, 1, 2, 5, 10, 30 secs and 1 min. Post-run manipulation includes rescaling, curve tracking and selection of the part of the curve required for the rate calculation. Rate is calculated using a linear regression algorithm before multiplying by the entered factor.






#### **Multi-Wavelength**

Up to 10 wavelengths may be entered, allowing the measurement of multiple wavelengths on a series of samples.



		Multi	-Wavelength Test Report
	Name: DFE.q and Time: 12	on 1.05.2008 10:0:	5:39
No	500,0pm	400.0om	300.0nm
		0.42	0.81



#### **Performance Validation**

for the GLP compliant laboratory SpectroQuest spectrophotometers may be automatically selfcalibrated on switch-on, using the 656.1 nm deuterium emission line. This function may be repeated at any time.

The wavelength accuracy may be checked using the "WL Validity" program (wavelength calibration standards required).

The absorbance accuracy at several wavelengths may be checked using the "Accu Validity" program.

Each SpectroQuest spectrophotometer includes all of the functionality shown on pages 20 and 21 as a standard feature. See pages 22 and 23 for optional software packages to expand your applications.




# Optional Advanced Windows-Based Application Software

The SpectroQuest Windows®-based PC application software takes the best features of the integrated operating software plus more powerful data processing and expanded data collection and storage capability. It comes standard with SpectroQuest PC models (SQ-2802PC and SQ-2802PCS) and is optional on all other SpectroQuest models.

#### The PC application software offers:


- 1. Abs/%T/Conc Test
- 2. Quantitative (standard curve)
- 3. Kinetics
- 4. Multi-wavelength Test
- 5. Wavelength Scanning
- 6. DNA/Protein

#### **Multi-Wavelength**



Up to 32 wavelengths can be selected and multiple samples can be measured. (Auto cell changer is required to run multiple samples automatically).

#### **Quantitative (Standard Curve)**




Use up to 32 standards to establish standard curve.

Four methods for fitting a curve:

- 1. Linear fit
- 2. Linear through zero
- 3. Square fit
- 4. Segmented

#### **Wavelength Scanning**



Automatically record peaks and valleys. Eight channels can simultaneously store up to 8 curves. Post-run manipulation and processing includes:

- 1. Re-scaling axes, curve
- 2. Smoothing, combination, zooming, overlap...
- 3. 1st to 4th derivative




#### **SpectroQuest**<sup>TM</sup> **Software**



The SpectroQuest Windows-based PC application software takes the best features of the integrated operating software plus more powerful data processing and expanded data collection and storage capability. It comes standard with SpectroQuest PC models (SQ-2802PC and SQ-2802PCS) and is optional on all other SpectroQuest models.


#### **DNA/Protein**



Concentration and DNA purity are quickly and easily calculated: Absorbance ratios 260nm/280nm with optional subtracted absorbance at 320nm.

DNA Concentration =  $62.9 \times A260 - 36.0 \times A280$ Protein Concentration =  $1552 \times A260 - 757.3 \times A280$ Other wavelengths and factors may be entered.

#### **Kinetics (Abs vs. Time)**



The Kinetics mode may be used for time course scanning or reaction rate calculations. Abs. vs. time graphs is displayed on the screen in real time.

Wait time, measurement time and time intervals may be entered.

Post-run manipulation includes re-scaling, curve tracking and selection of the part of the curve required for the rate calculation. Rate is calculated using a linear regression algorithm before multiplying by the entered factor.

#### Data Download Software (Cat. No. 2800-405)

The basic data download software is designed to download data stored in the spectrophotometer memory to a PC in a text file format for easy exporting into a spreadsheet.

Item No.	Description
SQ2800-401	Advanced Application Software for PC's
	(for SQ-2800 series)
SQ2800-405	Basic Data Download Software
	(for SQ-2800 and SQ-2802 models)
SQ3802-401	Advanced Application Software for PC's
	(for SQ-3800 series)
SQ4802-401	Advanced Application Software for PC's
	(for SQ-4800 series)

Se habla español. Visítenos en el Internet: www.unicosci.com/espanol o llámenos por teléfono al 609-240-5507.



#### **Optional SpectroQuest Accessories**

#### Test Tube Holder (Cat. No. SQ2800-101P)

Test tube holder kit for 8-20mm diameter test tubes. Includes universal base, V-type tube holder. The maximum tube height is 100mm.



SQ2800-101P

#### Water-Jacketed Cell Holder (Cat No. SQ2800-105P)

Water-jacketed single 10mm cell holder kit including universal base and one water-jacketed cell holder. It maintains desired temperature by circulating constant-temperature water from water bath (water bath required and not included).



SQ2800-105P



#### 4-Cell 100mm Long Path Cell Holder (Cat. No. SQ2800-102P)

Rectangular long path cell holder kit for 4 cells up to 100mm pathlength.

SQ2800-102P



#### SQ2800-106P

#### Micro Cell Holder (Cat No. SQ2800-106P)

Measure a sample with volume of 100uL using micro cell holder. The x-y adjustable mechanism is used to align cell with optical beam for optimized results.

#### 4-Cell 50mm Long Path Cell Holder (Cat No. SQ2800-102-50)

Rectangular long path cell holder for 4 cells up to 50mm pathlength.



SQ2800-102P-50

#### Peltier Unit (Cat No. SQ2800-107P)

Peltier unit for continuous temperature control from 15° to 40° C. The x-y adjustable mechanism is used to align cell with optical beam for micro cell setup. The temperature display resolution is 0.1° C. The unit consists of a controller and a thermoelectrically controlled cell holder and SQ panel.



SQ2800-107P



#### Cylindrical Cell Holder (Cat No. SQ2800-104P)

Cylindrical cell holder kit for single cell up to 100mm pathlength (20mm dia.). Includes universal base and one holder.

SQ2800-104P



#### **SpectroQuest**<sup>TM</sup> Accessories



SQ2800-108P

setup.

#### Ambient Sipper Unit (Cat No. SQ2800-108P)

Sipper system for single cell flow thru. The x-y adjustable mechanism is used to align cell with optical beam for micro cell flow cell setup. The sipper unit consists of a flow-thru controller with peristaltic pump and flow-thru front panel (flow cell and tubing not included).



SQ2800-121

#### 8-Position Auto Cell Changer (Cat. No. SQ2800-121)

Eight-position automatic cell changer designed for SQ-2800/2802/3802 series spectrophotometers.

#### Peltier/Sipper System (Cat. No. SQ2800-109P)

Peltier/Sipper system for a single cell flow thru and continuous temperature control from 15° to 40° C. The x-y adjustable mechanism is used to align cell with optical beam for micro flow cell setup. The temperature display resolution is 0.1° C. The unit consists of a Peltier/Sipper controller with peristaltic pump and a thermoelectrically temperature controlled cell holder with panel. The unit can be used a flow thru only or temperature control only.



SQ2800-109P

#### 6-position Auto Cell Changer (Cat. No. SQ4802-120)

Six-position automatic rotating cell changer designed for SQ-4802 series spectrophotometer.



SQ4802-120

Note: requires flow cell and proper tubing to complete flow thru setup.

Note: Requires flow cell and proper tubing to complete flow-thru

We inventory an extensive line of additional parts and accessories not shown here. If you do not see what you want or need, please contact us for more information.



SQ2800-122

#### Reflectance Measurement Attachment (5° incident angle) (Cat. No. SQ2800-122)

The technique of reflectance measurement is used for evaluation of materials relative to a reflectance surface. The minimum sample is (L) 30 x (W) 30 mm.



#### **UNICO Cuvets**


#### **Cuvets and Glassware for UNICO Spectrophotometers**

A large variety of cuvets are available including









S90-311Q

standard round, or square glass cells, quartz cells,	
long path cells, semi-micro and micro cells, short	
path cells, cylindrical cells and flow thru cells.	
•	

Item #	Previous Item #	Materials	Shape	Pathlength	Interior Width	Height	Capacity	Wavelength Range	Qty per Pack
Regular or Macro	Cuvets (10 mm p	athlength)							
S-90-304-1G	2100-304-1G	Glass	Square	10 mm	10 mm	45 mm	3.5 ml	335 to 900 nm	1
S-90-304G	2100-304G	Glass	Square	10 mm	10 mm	45 mm	3.5 ml	335 to 900 nm	2
S-90-309-1Q	2100-309-1Q	Quartz	Square	10 mm	10 mm	45 mm	3.5 ml	200 to 1200 nm	1
S-90-309Q	2100-309Q	Quartz	Square	10 mm	10 mm	45 mm	3.5 ml	200 to 1200 nm	2
S-90-301	2100-301	Glass	Tube	10 mm	10 mm	100 mm	4.0 ml	335 to 1100 nm	12
S-90-302P-100	2100-302P	Polystyr.	Square	10 mm	10 mm	45 mm	3.5 ml	335 to 1100 nm	100
S-90-302P-500	2100-303P	Polystyr.	Square	10 mm	10 mm	45 mm	3.5 ml	335 to 1100 nm	500
S-90-305P	2100-305P	Polystyr.	Square	10 mm	4 mm	45 mm	1.5 ml	335 to 1100 nm	500
Non-10 mm Path	length Cuvets								
S-90-308-50G	2100-308-50G	Glass	Cylindr.	50 mm	n/a	22 mm	14.1 ml	335 to 2500 nm	1
S-90-308-100G	2100-308-100G	Glass	Cylindr.	100 mm	n/a	22 mm	28.2 ml	335 to 2500 nm	1
S-90-310Q	2100-310Q	Quartz	Rectan.	5 mm	w/ stopper	48 mm	1.7 ml	170 nm-2500 nm	1
S-90-311Q	2100-311Q	Quartz	Rectan.	50 mm	w/ stopper	48 mm	17.5ml	170 nm-2500 nm	1
S-90-312Q	2100-312Q	Quartz	Rectan.	30 mm	w/ stopper	48 mm	10.5ml	170 nm-2500 nm	1
5-90-313Q	2100-313Q	Quartz	Rectan.	20 mm	w/ stopper	48 mm	7.0 ml	170 nm-2500 nm	1
S-90-314Q	2100-314Q	Quartz	Rectan.	100 mm	w/ stopper	48 mm	3.0 ml	170 nm-2500 nm	1
S-90-320G	2100-320G	Glass	Rectan.	20 mm	open top	45 mm	7.0 ml	335 to 2500 nm	1
S-90-321G	2100-321G	Glass	Rectan.	50 mm	open top	45 mm	17.5 ml	335 to 2500 nm	1
5-90-322G	2100-322G	Glass	Rectan.	100 mm	open top	45 mm	35.0 ml	335 to 2500 nm	1
5-90-330G	2100-330G	Glass	Rectan.	1 mm	10 mm	45 mm	0.35 ml	335 - 1000 nm	1
5-90-331Q	2100-331Q	Quartz	Rectan.	1 mm	10 mm	45 mm	0.35 ml	170 - 2500 nm	1
5-90-332G	2100-332G	Glass	Rectan.	2 mm	10 mm	45 mm	0.70 ml	335 - 2500 nm	1
S-90-333Q	2100-333Q	Quartz	Rectan.	2 mm	10 mm	45 mm	0.70 ml	170 - 2500 nm	1
S-90-335Q	2100-335Q	Quartz	Rectan.	3 mm	10 mm	45 mm	1 ml	170 - 2500 nm	1
Flowcell, Macro	and Micro								
S-90-340FG	2100-340FG	Glass	Square	Frosted	10x7x40 mm	48 mm	3.0 ml	320 to 2500 nm	1
S-90-341FG	2100-341FG	Glass	Square	Frosted	10x4x40 mm	48 mm	1.8 ml	320 to 2500 nm	1
S-90-342FQ	2100-342FQ	Quartz	Square	Frosted	10x7x40 mm	48 mm	3.0 ml	170 to 2500 nm	1
S-90-343FQ	2100-343Q	Quartz	Square	Frosted	10x4x40 mm	48 mm	1.8 ml	170 to 2500 nm	1
S-90-344FG	2100-344FG	Glass	Square	Black	10x4x12 mm	48 mm	0.48 ml	320 to 2600 nm	1
S-90-345FG	2100-345FG	Glass	Square	Black	10x3 mm (dia.)	48 mm	0.07 ml	320 to 2500 nm	1
S-90-346FQ	2100-346FQ	Quartz	Square	Black	10x4x12 mm	48 mm	0.48 ml	170 to 2500 nm	1
S-90-347FQ	2100-347FQ	Quartz	Square	Black	10x3 mm (dia.)	48 mm	0.07 ml	170 to 2500 nm	1
Semi-Micro, Sub-	Micro Cuvets								
S-90-350G	2100-350G	Glass	Lid	Black	10x4	45 mm	1.0 ml	320 to 2,600 nm	1
S-90-351Q	2100-351Q	Quartz	Lid	Black	10x4	45 mm	1.0 ml	170 to 2,600 nm	1
S-90-353G	2100-353G	Glass	Lid	Black	10x2	45 mm	0.5 ml	320 to 2,600 nm	1
5-90-354Q	2100-354Q	Quartz	Lid	Black	10x2	45 mm	0.5 ml	170 to 2,600 nm	1
S-90-356G	2100-356G	Glass	Lid	Black	10x1	45 mm	0.25 ml	320 to 2,600 nm	1
S-90-357Q	2100-3 <i>57</i> Q	Quartz	Lid	Black	10x1	45 mm	0.25 ml	170 to 2,600 nm	1
5-90-358Q	2100-358Q	Quartz	Stopper	Black	10x2 x 5 mm	45 mm	0.100 ml	170 to 2,600 nm	1
S-90-359Q	2100-359Q	Quartz	Stopper	Black	10x2 x2.5 mm	45 mm	0.050 ml	170 to 2,600 nm	1
3-70-33 <b>7</b> 4	Z100-334Q	Quariz	Siopper	DIGCK	TUXZ XZ.J MM	45 mm	0.030 mi	170 10 2,000 nm	I



#### **UNICO Laboratory Products**

# Please see these and additional laboratory products from UNICO on our website at www.unicosci.com.



PowerSpin™ C822 Model HX Centrifuge



PowerSpin™ C8606 Model MX Horizontal Centrifuge



PowerSpin™ C8624 Model MX Centrifuge,



MTR22 Programmable Multi-Mixer



LVM 1000 Vortex Mixer



LTTR200 Rock-IT™ Tube Mixer



G388 Dual Binocular Head Infinity Optics Microscope



**IP750 Infinity Microscope** 



ZM180 Zoom Stereomicroscope



**Glove Box Dispenser** 



Tube-CUBE™ Tube Rack



Phlebotomy Tray



# 











#### **Contact Us:**

UNICO
United Products and Instruments, Inc.
182 Ridge Road, Suite E
Dayton, New Jersey 08810
U.S.A.

Toll Free Phone: 800-588-9776 Direct Dialing: 732-274-1155 Fax: 732-274-1151

E-mail: info@unicosci.com Website: www.unicosci.com

En español: 609-240-5507 E-mail: ventas@unicosci.com