



# PowerSpin™ MX Horizontal Centrifuge

UNICO® is proud to introduce our powerful 6-place Horizontal Spin centrifuge, the new PowerSpin™ MX! The PowerSpin™ MX spins 1000 – 3400 RPM, to insure a clean and complete separation of cellular components yielding superior serum or plasma specimen. Horizontal swing-out rotor provides high G-force and flat separation meniscus line.

The all digital push button front panel makes it easy to select or modify speed and run time even while wearing gloves. The PowerSpin™ MX also features a two line digital display that shows status, time, speed (RPM), and RCF (g).

Each PowerSpin™ Series centrifuge runs extremely quiet. The rubber shock absorbing suction cup feet hold the unit in place resulting in a quieter run with less vibration

Integrated True Zero RPM Safety System prevents lid from being open while the rotor is in motion.

Call us today for more information.

Tube Holders are  
“Vertical at Rest”  
to allow use of  
Cover Slips for  
Fecal Protocols



### Features and Specifications:

- 6 place swing-out horizontal spin rotor for 2 to 15 ml tubes
- Push button panel for convenience when wearing gloves
- Quiet running , only 50 decibels at 3400 RPM
- Stainless steel rotor chamber is accessible for easy cleaning
- Easy to read digital display shows real time RPM and RCF readings
- Remembers last setting for quick and easy repeat testing
- Variable Speed, 1000-3400 RPM
- Suction cup rubber feet eliminate vibration
- 60 minute digital timer with tone signal
- On/off switch located on rear panel of unit
- Detachable power cord
- Lid window for speed verification with strobe tachometer
- Zero RPM Safety Latching System to prevent lid from being opened while rotor is spinning
- Manual lock override in case of power outage
- Indicator light signals when the power is on
- Standard 115V/60Hz AC (220V/50Hz available on request)

### Ordering Information:

Model	Catalog #	Speed Range	Speed Monitor	Max. RCF	Capacity	Swing-out Rotor	Motor
<b>PowerSpin™ MX</b>	<b>C8606</b>	1000 - 3400 RPM	Digital display	1825 G	6 x 2-15 ml	6 place	AC 115V

### Accessories and parts:

C8606-01	6 place swing out rotor	Centrifuge Size:	13 ¾ x 15 x 11”
C800-02	Regular tube shield, pack of 2	Shipping Size:	18 in x 19 in x 15.5”
C800-03	Adapters for 5 ml tubes, pack of 10	Net Weight:	31 lbs
C800-04	Pediatric or 2 ml tube adapters, pack of 2	Shipping Weight:	36 lbs
C800-10	Fuse 250A 4A, each	Warranty:	One year, parts and labor