




The top shelf of each T-Rack has cross-referencing Alpha – Numeric designations so you can track the exact location of specific tubes in the rack.



Four T-Racks shown joined together, both 10-13mm and 13-16mm racks are attached. Mix and match to meet your needs, or to help with Accessioning.



Mix and match to meet your needs, or to help with Accessioning.

The only rack you will ever need! Starting with a single 50 place T-Rack, you can attach as many racks together as you have room for! These patented racks quickly and easily attach together with their unique built-in T-Connector, which allows you the flexibility to attach one, two, three, or however many racks you decide together!

How do I attach T-Racks together?

It's so easy, even a hospital administrator can do it!

T-Rack's T-Connectors have a positive "click-stop", which locks the racks together, so they will not slide apart while being handled, or moved around the bench top, or refrigerator, etc. This click-stop holds them tightly together, yet is very easy to disassemble as well!

To disassemble, simply slide the racks apart and use as single rack, or join together with other T-Racks.

You can attach 10-13 mm racks to 13-16 mm racks! Attach any configuration to fit your needs! T-Racks allow you the flexibility to select how best to use them in your lab.

T-Racks can be used to help with Accessioning and are a great batching aide.

Item No.	Description	UOM
54500	T-Rack 50 Place, 13-16 mm, Expandable	Each
54500-2	T-Rack 50 Place, 13-16 mm, Expandable	Pk/2
54510	T-Rack 50 Place, 10-13 mm, Expandable	Each
54510-2	T-Rack 50 Place, 10-13 mm, Expandable	Pk/2
54520	T-Rack 50 Place, (25 holes of 13-16 mm & 25 of 10-13 mm), Expandable	Each
54520-2	T-Rack 50 Places, (25 holes of 13-16 mm & 25 of 10-13 mm), Expandable	Pk/2

Specifications:

- Each T-Rack holds up to 50 tubes in two 25 hole sections
- The top shelf of each T-Rack has cross-referencing Alpha – Numeric designations so you can track the exact location of specific tubes in the rack
- Attach additional T-Racks together, make a single large rack, as big as you need!
- Made of durable, autoclavable polypropylene
- T-Racks will not break, or become brittle in low temperature refrigerators
- Easy to clean and disinfect
- Dimensions: Each 50 Place T-Rack measures: 8.75" x 3.75" x 5.25"
- Your T-Racks will arrive disassembled and snap them together in seconds

Storage Tray



Organize and store supplies, tube racks (as shown) and more in this three-section, see-thru tray. Great for batching samples and organizing your work flow, or to store materials in the refrigerator, etc. Made of durable 3/16" acrylic.

Item No.	Description	UOM	Size
52030-1	Storage Tray, 3 Section	EA	12.75"W x 10"D x 3"H