

- Twenty-four (24) place rotor with replaceable rubber ring and screw down metal lid.
- Nominal speed is 12,000 RPM
- Dependable mechanical 5-minute timer with bell signal at end of run
- Micro-Hematocrit reader included
- Air-cooled, brush type AC motor
- All metal housing for years of dependable use
- Safety switch in latch shuts off motor when lid is lifted
- Electric motor braking
- Indicator light tells when unit is in use



UNICO Part #	Description	UOM
C-MH30	24 place Micro-Hematocrit Centrifuge with Reader	each
C-MH30-Reader	Micro-Hematocrit Reader	each

UNITED PRODUCTS & INSTRUMENTS

182 Ridge Road, Suite E Dayton, New Jersey 08810 * Tel: 800-588-9776 * Fax: 732-274-1151 * E-mail: info@unicosci.com * Web: www.unicosci.com

© 2011 UNICO®