

# UNICO® Colpo-Master™ I Colposcopes

## Zoom Optics with LED, or Halogen Illumination

UNICO® Colpo-Master™ I Suspension Arm Colposcopes feature continuous zoom 4 – 27x magnification with 20x High Eyepoint eyepieces and sharp focus for 3-D imaging. Available with Binocular or video-ready Trinocular heads. Video-ready heads accept C-mount cameras; we offer both Digital and Analog S-Video options, EMR compatible.

Designed with your viewing comfort in mind; each model has dual diopter adjustments to independently adjust the focus of both Left and Right eyepieces to eliminate eyestrain; all with 298mm working distance. The field of view (FOV) ranges from 58.5mm at 4x to 8.5mm at 27x. All models with your choice of 3 or 5 leg sturdy mobile bases.

Your choice of LED or Halogen illuminated models. Super bright LED models feature dual rheostats, one each for white and green LED intensity control.

Colpo-Master™ I's Suspension Arm is rock-steady. To make head positioning easy it comes with locking casters and stainless steel handlebars, standard.

Complete with dust cover, power cord, manual, tool kit and 7 Year Limited Warranty, 1 year on electrical system (excludes bulbs and fuses). Halogen models include a spare bulb. Optional accessories are available - please contact UNICO for additional information.



- Featuring 4 White and 1 Green LED's
- Both White and Green LED offering individual brightness adjustments
- 30,000 hours LED Bulb Life

CS-103LED

UNICO ITEM #	DESCRIPTION COLPO-MASTER™ I	U/M
CS-103	Colpo-Master™ I Binocular Colposcope, Halogen, 3 Leg	1 ea
CS-103LED	Colpo-Master™ I Binocular Colposcope, LED, 3 Leg	1 ea
CS-103T	Colpo-Master™ I Trinocular Colposcope, Halogen, 3 Leg	1 ea
CS-103T-LED	Colpo-Master™ I Trinocular Colposcope, LED, 3 Leg	1 ea
CS-105	Colpo-Master™ I Binocular Colposcope, Halogen, 5 Leg	1 ea
CS-105LED	Colpo-Master™ I Binocular Colposcope, LED, 5 Leg	1 ea
CS-105T	Colpo-Master™ I Trinocular Colposcope, Halogen, 5 Leg	1 ea
CS-105T-LED	Colpo-Master™ I Trinocular Colposcope, LED, 5 Leg	1 ea
CS-9001	Halogen Replacement Bulb, 150 W	1 ea

United Products & Instruments, Inc  
 182 Ridge Road, Suite E Dayton NJ 08810  
 800-588-9776 \* [info@unicosci.com](mailto:info@unicosci.com) \* [www.unicosci.com](http://www.unicosci.com)

Colpo-Master™ is a trademark of United Products & Instruments, Inc, dba UNICO

