

# DermaScout™ Dermatoscope



- 6x Magnification
- 9 LED Illumination
- Measuring Scale (mm)

- The DermaScout is efficient and dependable instrument designed for busy professionals
- Helpful in distinguishing benign from malignant lesions, early detection of melanoma
- View and measure skin lesions, fungal infections, skin tumors, warts and other conditions



DermaScout™ shown in storage case

## Features:

6x Magnification, with 40 mm top lens

40 mm Field of View, (flat field)

LED Illuminated with an array of 9 LED's

Reticule, Dual Measuring scales: mm & inches

Measures in 0.1 mm increments

Storage / Transport Case

Diopter focus control for sharp focusing

31 mm Working Distance

Battery powered, over 150 hours

20 mm: 10 left of center/10 right of center

Measures in .005" increments

Battery included (Type CR2, 3 Volt)

Item #	Description	UOM
DS-6	DS-6 DermaScout™ dermatological skin magnifier	Each

# UNICO®

[www.unicosci.com](http://www.unicosci.com) \* [info@unicosci.com](mailto:info@unicosci.com)

Toll Free: (800) 588-9776 \* Phone: (732) 274-1155 \* Fax: (732) 274-1151