

# PowerSpin™ EX Centrifuge

UNICO® is proud to present the PowerSpin™ EX Centrifuge. With 24 place horizontal rotor this centrifuge is a great answer for a busy laboratory.

The PowerSpin™ DX has variable speed 500 – 3500 RPM, to insure a clean and complete separation. Quiet and dependable, with a maintenance free motor, with stainless steel tube shields it can be used to process blood, urine, fecal and other specimens.


The all digital front panel with 3 pre-programmed options

- 3200 RPM for 10 minutes for Blood
- 1800 RPM for 5 minutes for Urine
- 1300 RPM for 10 minutes for Fecal

Custom setting are easy and intuitive with all digital panel and oversized display.



Model #C8924



- Zero RPM Safety Latching System prevents lid from being opened while rotor is spinning
- Quiet running, only 50 decibels at 3500 RPM
- 99 minute digital timer, 1 sec. increment
- 500-3500 RPM, 100 RPM increment

- 3 pre-programmed options for blood, fecal and urine
- Remembers last setting for quick repeat runs
- Adjustable Acceleration and Deceleration
- Stainless steel rotor chamber for easy cleaning
- Available in 110V and 220V - Contact UNICO for details

MFG Part #	Description	UOM
C8924	PowerSpin EX, 24 place horizontal spin centrifuge. 110V	EA
C8924E	PowerSpin EX, 24 place horizontal spin centrifuge. 220V	EA