

FREEDOM JR.™

Freedom Jr.™ Microscope Features

Free from cords

- Standard AA batteries

Free from charging downtime

- Over 240 Hours with a single battery change

Free from Lamp Changes

- New LED design gives a lamp life over 100,000 hours

Free of worries

- Student-proof design
- Free from the risk of A/C current
- Reverse nosepiece
- All metal construction
- DIN Standard - High grade optics
- Free to use on any work surface
- Free from heat
- Limited lifetime warranty

New Illumination System New user friendly conveniences



You told us what you wanted and **UNIC** has delivered. A portable microscope that has brilliant illumination and safe from risky high voltage & A/C current. Our answer: **UNIC FREEDOM Jr.™**
CORDLESS MICROSCOPE!

Specifications

- Reversed positioned triple nosepiece with click stop.
- Widefield 10x eyepiece, rotatable with safety lock and built-in pointer.
- Monocular head inclined at 45 degrees.
- Large stage, acid and reagent, resistant with lock-on stage clips and safety stop.
- Stage drilled to accept optional attachable mechanical stage.
- Intermediate focus adjustment (eliminates the complexity of making two adjustments).

- Built-in N.A. 0.65 condenser with 5-position disc diaphragm.
- Optional A/C Adapter
- Plastic dust cover and user's manual.

Actual Size : 6 1/2" x 4 1/2" x 11 1/4"
Shipping Dimension: 8" x 5 3/4" x 13"
Shipping weight: 4.5 lbs.



United Products & Instruments, Inc.

162 Ridge Road, Suite E, Dayton, NJ 08810 USA Tel: 800-588-9776 / 732-274-1155 Fax: 732-274-1151
E-MAIL: sales@unico1.com Visit us at www.unico1.com