

# UNICO® Count On - It!

## 2, 5 and 8 key Differential Counters



Manual differential counters have been used for decades in the clinical laboratory. Models L-BC9 and L-BC6 featuring a WBC maturation series picto-strip just above the key windows. **Unico®** offers the three most popular counters: 2 key, 5 key with total window, and 8 key with total window

### Features:

- Three models: 2 key, 5 key and 8 key
- WBC maturation series Picto-strip
- Heavy duty, impact resistant plastic housing
- Dual reset knobs
- Wide base with rubber feet for stability
- Signal bell when 100 cells have been counted

	L-BC3 – 2 key	L-BC6 – 5 key	L-BC9 – 8 key
Size	4 in x 3.0 in x 2.0 in	9.7 in x 3.0 in x 2.0 in	12.3 in x 3.0 in x 2.0 in
Shipping Size	5.5 in x 4.5 in x 3.7 in	11.3 in x 4.5 in x 3.7 in	13.9 in x 4.6 in x 3.6 in
Weight	0.2 lbs	1.0 lbs	1.4 lbs
Shipping Weight	0.5 lbs	1.4 lbs	1.8 lbs
Limited Warranty	One year	One year	One year