


UNICO®

Shake-It!

6-position Pipette Shaker



This type of shaker has been a reliable standard in the laboratory for generations. The Shake-It is ideal for use in hematology, for special tests and manual techniques.

Features

- Room for 6 pipettes
- 15 min. timer
- Triple position switch: On/Off/Continuous use
- Power indicator light
- All metal construction
- Acid/ reagent resistant paint finish
- 110V AC

Size	5.2 in x 5.8 in x 6.2 in (133 mm x 147 mm x 157 mm)
Shipping Size	8.1 in x 6.9 in x 8.1 in (206 mm x 176 mm x 206 mm)
Weight	3.8 lbs (1.7 kg)
Shipping Weight	4.6 lbs (2.1 kg)
Limited Warranty	One year

United Products & Instruments, Inc. 182 Ridge Road, Suite E, Dayton, New Jersey, 08810, USA
Tel: 1-800-588-9776 Fax: 732-274-1151 E-mail: sales@unicol.com Web Site: <http://www.unicol.com>