

# UNICO®

# Rotate-it!



## Variable Speed Laboratory Rotator

The UNICO® L-RT28 rotator is a compact, durable, variable speed orbital mixer. Designed to accommodate a wide variety of printed or glass well slides, culture plates, flasks and small test tube racks, as well as RPR and VDRL slides or 96 well microtiter (immunoassay) plates.

### Features

- 11" x 11" mixing plate with removable rubber pad
- Variable speed 25-250 rpm
- 30 min. timer
- Orbital stroke of 20 mm
- 3 position switch: Constant on/off/timed mixing
- Analog tachometer (adjustable)
- Small flask hold down spring included
- Holding one 1000 ml flask, two 500 ml flasks, and four 250 ml flasks
- Limited One year Warranty

Size	10.4" x 9.2" x 6.0"
Shipping Size	12.2" x 11.2" x 7.4"
Weight	9.8 lbs
Shipping Weight	10.8 lbs

### **United Products & Instruments, Inc.**

182 Ridge Road, Suite E, Dayton, New Jersey, 08810, USA  
Tel: 1-800-588-9776 \* Fax: 732-274-1151

E-mail: [sales@unicol.com](mailto:sales@unicol.com) \* Web Site: [www.unicol.com](http://www.unicol.com)