



Rock-IT!



Rock-IT Jr!

The **UNICO®** models **LTTR 100** and **LTTR 200 Rock-IT** Tube Rockers are heavy-duty, yet nearly silent mixers designed for years of everyday use. Smooth and gentle movement at 25 RPM assures a uniform, thoroughly mixed specimen, while preserving the integrity of the fragile blood cells. Rock-IT's will also get hematology controls ready to use in no time for accurate results. The reversible non-slip pad will accommodate larger and hard to mix tubes. Rock-IT's are an economical choice for your lab; let the Rock IT work for you!

- Capacity of 8 or 15 tubes, up to 120 mm long
- 25 rocking cycles per minute, for gentle mixing
- Non-slip, reversible and autoclavable tube pad, silicon rubber
- Compact size uses minimal benchtop space
- One year limited warranty

Item Number	Description	
LTTR100	Rock-IT Jr. Tube Rocker,	8 Place
LTTR200	Rock-IT Tube Rocker,	15 Place
LTTR100-1	Replacement Tube Pad,	8 Place
LTTR200-1	Replacement Tube Pad,	15 Place
Size and Shipping Info	LTTR 100	LTTR 200
Size	8.6" x 4.7" x 5.0"	13.7" x 4.7" x 5.0"
Shipping Size:	10.6" x 6.7" x 6.9"	15.7" x 6.7" x 6.9"
Net Weight / Shipping Wt.:	4.0 lbs / 5.0 lbs.	4.8 lbs / 5.8 lbs.