

M1 and M2 Prism Microscopes

This amazing monocular microscope is truly portable and economical! Simple and easy to use, students adapt quickly to this microscope and enjoy discovery. No electricity required; the unique prism gathers and directs ambient room light and brightly illuminates slides and other specimens.

Focus is easily controlled by simply sliding the monocular tube up and down in its holder and the microscope can pivot on its base allowing the user to incline the viewing tube for viewing comfort or to gather more light.

Model M1: The combination of a 4X objective and 10X eyepiece provide 40X power

Model M2: The combination of a 4X objective and 15X eyepiece provide 60X power

Includes a locked-on eyepiece and stage clips, a dust cover and instruction manual. Limited Lifetime Warranty.



Item No.	Description
M1	Compound Microscope, 40x
M2	Compound Microscope, 60x

Actual Size: 4.75" x 4" x 10"

Shipping Dimension: 6" x 5.25" x 11.5"

Net Weight: 2 lbs

Shipping Weight: 2.5 lbs

	M1	M2
Head	Monocular	Monocular
Eyepiece	10X WF	15X WF
Objectives	4X	4X
Condenser	None	None
Diaphragm	None	None
Stage	Plain	Plain
Illumination	Prism	Prism
Bulb	None	None
Focusing	Tube	Tube

M1 Prism Microscope

United Products & Instruments, Inc.
800-588-9776 www.unico1.com

UNICO[®]