

High School Microscopes

M230 FREEDOM Rechargeable Series


Built on the reliable and time tested M220 frame, the M230 series incorporates all the great features of the M220's with the added versatility of a rechargeable battery operated system and an adjustable rheostat light intensity control. Use this model as either a portable field-use scope on battery power, or plugged into a wall outlet for traditional classroom use with the AC adapter/recharger. Recharging system uses 3 rechargeable AA batteries (included).

Available with or without add-on graduated mechanical stage. Each microscope includes AC adapter/recharger, dust cover and instruction manual. Limited Lifetime Warranty: Excludes bulbs and batteries.

Cordless LED Rechargeable Microscopes

Item No.	Description
M230	High School Monocular Microscope, LED Illumination, Rechargeable
M230-M	High School Monocular Microscope, LED Illumination, Rechargeable, Mechanical Stage
B2-5101	Graduated Mechanical Stage, Attachable

	M230	M230-M	M230-D	M230-DM
Head	Monocular	Monocular	Dual head	Dual Head
Eyepiece	10X WF	10X WF	10X WF (2)	10X WF (2)
Objectives	4X, 10X, 40X	4X, 10X, 40X	4X, 10, 40X	4X, 10, 40X
Condenser	NA 0.65	NA 0.65	NA 0.65	NA 0.65
Diaphragm	Disc, 5-hole	Disc, 5-hole	Disc, 5-hole	Disc, 5-hole
Stage	Plain, tapped	Mechanical	Plain, tapped	Mechanical
Illumination	Battery/AC	Battery/AC	Battery/AC	Battery/AC
Bulb	LED	LED	LED	LED
Focusing	Coarse/Fine	Coarse/Fine	Coarse/Fine	Coarse/Fine



M230 Rechargeable Microscope

Actual Size: 8.5" x 5" x 11.5"
Shipping Dimension: 12" x 8" x 16"
Shipping Weight: 9 lbs