

# ST60 SERIES

## Stereo Microscope

This dual magnification stereo microscope features high-resolution optics and precise optical alignment that assures crisp sharp and flat 3- dimensional image. UNICO cool lite illumination system gives you optimum light at ambient temperature to prolong life of live material. It is an ideal instrument for small dissection, crystal analysis, quality inspection and all educational applications.

### Specifications

- Paired, widefield 10X/FN 18mm eyepieces
- Interpupillary distance adjustable from 51mm to 75mm
- Diopter adjustment on right eyetube for maximum acuity
- Inclined binocular body with choice of magnifications: 1X/2X, 1X/3X or 2X/4X
- Binocular body rotates 360 degrees
- 65mm working distance
- Smooth accurate focusing by rack and pinion with tension control and slip clutch
- Built-in dual illumination with bright white fluorescent bulbs
- Four position switch: off, reflected light, transmitted light, or reflected and transmitted light simultaneously
- 75mm diameter frosted stage plate
- Plastic cover
- 5 year limited warranty



MODEL	EYEPIECE	OBJECTIVE	ILLUMINATION	MAGNIFICATION
ST62FF	WF10X	1X/2X	Fluorescent	10X/20X
ST63FL	WF10X	1X/3X	Fluorescent	10X/30X
ST68FL	WF10X	2X/4X	Fluorescent	20X/40X

**United Products & Instruments, Inc.**  
*Manufacturers of Quality Optical and Analytical Systems*  
182 Ridge Road, Dayton, NJ 08810 USA  
Tel: 800-588-9776 ■ 732-274-1155 ■ Fax: 732-274-1151  
E-MAIL: sales@unico1.com  
WEB SITE: www.unico1.com

Shipping Dimensions: 370mm x 290mm x 230mm  
(14.6in. x 11.4 in. x 9 in)  
Shipping Weight: 10 lbs (5.6 kgs)

