


# Tube Adapter

Item: 8604-02



The 8604-02 Tube Adapter is designed for use with the C8704 centrifuge, enhancing its versatility by accommodating smaller tube sizes. This adapter is tailored for 5-15ml tubes with a length range of 100mm to 120mm, making it ideal for a variety of laboratory applications that require precise and efficient sample processing.

## COMPATIBILITY

<b>Centrifuge Model</b>	C8704
<b>Rotor Type</b>	4 Place Horizontal Rotor
<b>Tube Size</b>	5-15ml
<b>Tube Length Range</b>	100mm - 120mm