

PowerSpin™ HX Series Centrifuge

The horizontal spin fixed speed PowerSpin™ HX centrifuge spins at 3400 RPM. The horizontal rotor ensures a clean, complete separation of cellular components. Using a horizontal rotor for spinning Serum Separator Tubes (SSTs) creates a tight uniform horizontal seal, reducing the risk of red cell leakage. The HX series features a Zero RPM safety lock, so the lid can not be opened while the rotor is in motion. A digital timer and LCD display provide users with the precise control of the spin cycle. Powerful brushless motor is maintenance free for years of dependable use. Stainless steel guard bowls are easy to clean and maintain.



C822-01



C8300-03



Actual Size: 14.4" x 11.6" x 8.9"
Shipping Dim.: 16.5" x 13" x 10"
Net Weight: 22 lbs

Key Features:

- **Glove-Friendly Keypad Controls:** Easy-to-use push-button keypad designed for operation while wearing gloves.
- **Versatile Rotor Options:** Available with 6, 8, or 12-place rotors for standard 3 to 15 mL tubes. (Adapters may be required for tubes smaller than 10 mL or shorter than 100 mm.)
- **Variable Speed Control:** Adjustable speed range from 500 to 3500 RPM.
- **Programmable Timer:** Electronic timer adjustable from 1 to 99 minutes with an audible end-of-run signal. Automatic memory recall of the last run's time and speed settings simplifies repetitive protocols.
- **Memory Function:** Retains the last used time and speed settings even after the unit is powered off. These settings are automatically recalled upon restart.

C8312



C8306

Item No.	Description
C8302	HX, 6 Place Horizontal Spin Centrifuge
C8304	HXV, 12 Place Fixed Speed Centrifuge
C8306	HXV, 6 Place Variable Speed Horizontal Centrifuge
C8308	HXV, 8 Place Variable Sped Hoizontal Centrifuge
C8312	HXV, 12 Place Variable Speed Centrifuge
C8302-02	6 Place Swing Out Rotor
C8300-03	12 Place Angled Rotor
C8300-10	Fuse for HX Series Centrifuge
C822-02	Tube Shield, Pkg. of 2
C800-03	Short Tube Adapters, for 4-7 ml tubes, Pk/10
C800-04	Small Tube adapters, for 2 ml tubes, Pk/2
C800-34	Microhematocrit Adapters, 4 tubes each, Pk/2