



unicosci.com

182 Ridge Road, Suite E, Dayton, NJ 08810 USA

Email: info@unicosci.com

Phone: 732-274-1155

## UNICO® PowerSpin™ LX Series Centrifuge

**Maximum Speed:** 300-4000 RPM

**Maximum RCF:** 1860 xg

**Capacity:** 15ml × 8 tubes

**Noise Level:** ≤60 dB(A)

**Drive Motor:** Brushless DC

**Power Supply:** 100V-240V, 50Hz/60Hz

The UNICO PowerSpin LX Centrifuge combines high performance with versatile functionality, making it an ideal choice for various laboratory applications. Designed with a robust brushless DC motor and user-friendly controls, this centrifuge offers reliable and efficient separation for a range of sample types. To further streamline operations, the centrifuge features multiple preset programs specifically designed for common sample types, including **blood**, **urine**, and **fecal** samples. These pre-configured settings enhance accuracy and reliability, allowing users to achieve consistent results without the need for extensive manual adjustments.



**Model: LX-C85608**

LX Series Centrifuge is TÜV/UL certified and complies with UL, CAN/CSA, IEC safety, and FCC electromagnetic compatibility standards.

Item	Description
LX-C85608	LX, Variable Speed, 8 Place Metal Rotor Centrifuge
C8000-34	Microhematocrit Tube Adapter Set, 4 Place, Set of 2
C8000-03	Short Tube Adapter, 3-5ml (75mm) Blood Tubes, Set of 2