

UNICO® PowerSpin™ LX PRO Centrifuge

Maximum Speed: 300-4000 RPM

Maximum RCF: 1860 xg

Capacity: 15ml × 8 tubes

Noise Level: ≤60 dB(A)

Drive Motor: Brushless DC

Power Supply: 100V-240V, 50Hz/60Hz

The LX PRO Centrifuge combines high performance with versatile functionality, making it an ideal choice for various laboratory applications. Designed with a robust brushless DC motor and user-friendly controls, this centrifuge offers reliable and efficient separation for a range of sample types. It features preset programs for **blood**, **urine**, and **fecal** samples, enhancing convenience and accuracy in laboratory settings.

Model: LX-C85608



C800-02-1



Regular Tube shields for FX/LX/DX centrifuges. 1 piece.

Certified by TÜV Rheinland, this centrifuge meets rigorous safety and EMC standards for reliable laboratory use.



UNICO® PowerSpin™ FX PRO Centrifuge

Maximum Speed: 3500 rpm

Maximum RCF: 1397×g

Capacity: 15ml x 8 tubes

Noise Level: ≤60 dB(A)

Drive Motor: Brushless DC

Power Supply: 110V, 60Hz, 2A

The FX PRO Centrifuge is a high-performance laboratory instrument designed to deliver accurate and reliable results for a wide range of applications. With its advanced technology and user-friendly features, the FX Centrifuge provides a perfect balance of speed, capacity, and safety. Ideal for clinical, research, and quality control environments, this centrifuge ensures efficient sample processing and optimal results in a compact, durable package.

Model: FX-C80508



C8000-02



15ml Adapter, 24 pack available.

Certified by TÜV Rheinland, this centrifuge meets rigorous safety and EMC standards for reliable laboratory use.



UNICO® PowerSpin™ PX PRO Centrifuge


Maximum Speed: 3500 rpm
Maximum RCF: 1397×g
Capacity: 15ml x 8 tubes
Noise Level: ≤60 dB(A)
Drive Motor: Brushless DC
Power Supply: DC 12V, 3A

Elevate your laboratory's efficiency and precision with the UNICO PX Centrifuge, a state-of-the-art instrument designed to meet the demands of diverse scientific and clinical applications. Combining advanced technology with user-friendly features, the UNICO PX delivers exceptional performance.

Model: PX-C82608



C800-06



Short tube adapter for 3-5ml (75mm) blood tubes. Set of 2 pieces.

Certified by TÜV Rheinland, this centrifuge meets rigorous safety and EMC standards for reliable laboratory use.



UNICO® PowerSpin™ PXV PRO Centrifuge

Maximum Speed: 500 to 3800 rpm
Maximum RCF: 1647×g
Capacity: 15ml x 8 tubes
Noise Level: ≤60 dB(A)
Drive Motor: Brushless DC
Power Supply: DC 12V, 3A

Experience cutting-edge centrifugation technology with the UNICO PXV Centrifuge, a high-performance laboratory instrument designed to deliver superior results across a variety of scientific and clinical applications. The PXV combines robust functionality with advanced features, making it an essential tool for precise sample processing in any laboratory setting.

Model: PXV-C82708



C800-04



Adapter for pediatric or microtube (38-45mm tube length), pack of 2

Certified by TÜV Rheinland, this centrifuge meets rigorous safety and EMC standards for reliable laboratory use.

