



182 Ridge Road, Suite E, Dayton, NJ 08810 USA

Email: info@unicosci.com

Phone: 732-274-1155

unicosci.com

## UNICO® PowerSpin™ PX Centrifuge

**Maximum Speed:** 3500 rpm

**Maximum RCF:** 1397×g

**Capacity:** 15ml x 8 tubes

**Noise Level:** ≤60 dB(A)

**Drive Motor:** Brushless DC

**Power Supply:** DC 12V, 3A

Elevate your laboratory's efficiency and precision with the UNICO PX Centrifuge, a state-of-the-art instrument designed to meet the demands of diverse scientific and clinical applications. Combining advanced technology with user-friendly features, the UNICO PX delivers exceptional performance.

The PowerSpin™ PX Centrifuge comes standard with a Vehicle Power Adapter Cable.



**Model: PX-C82608**

The **C826-12** AC Power Adapter unit with a UL power cord is available for purchase.

This centrifuge meets top safety and electromagnetic compatibility standards, including **UL 61010-1**, **IEC 61010-2-020**, and **FCC Part 15**.



**C826-12**