

182 Ridge Road, Suite E, Dayton, NJ 08810 USA

Email: info@unicosci.com Phone: 732-274-1155



182 Ridge Road, Suite E, Dayton, NJ 08810 USA

Email: info@unicosci.com
Phone: 732-274-1155

UNICO® PowerSpin™ DX PLUS Series Centrifuge

Maximum Speed: 4000 RPM

Maximum RCF: 2508×g Capacity: 6 x 15ml tubes Noise Level: ≤62 dB(A) Drive Motor: Brushless DC

Power Supply: 100V-240V, 50/60Hz

Dimensions: 364 (L) × 440 (D) × 268 (H) mm

Weight: 14.5 kg

The UNICO DX Series Centrifuge stands out as an essential tool for precise and reliable centrifugation in laboratory settings. Engineered for versatility and safety, this centrifuge offers a maximum rotational speed of 4000 RPM and an exceptional relative centrifugal acceleration of 2508×g, catering to a range of separation needs from blood to fecal samples.

Model: DX-C87504 DX-C87506 DX-C87524

The PowerSpin DX PLUS comes with 3 Rotor options:

4 Place for 50ml tubes

6 Place Horizontal for 2-15ml tubes

24 Place Fixed Angle for 2-15ml tubes



4 Place Horizontal



6 Place Horizontal



24 Place Fixed Angle



UNICO® PowerSpin™ HX Series Centrifuge

Maximum Speed: 3400 RPM Maximum RCF: 1777 xg

Capacity: 15ml × 6 tubes
Noise Level: ≤ 62 dB(A)
Drive Motor: Brushless DC

Power Supply: 100V-240V, 50Hz/60Hz Dimensions: 364 (L) × 440 (D) × 268 (H)

Weight: 14.5 kg

Elevate your laboratory's efficiency with the UNICO HX Series
Centrifuge, designed for precision and reliability. Perfect for clinical and research applications, this high-performance centrifuge delivers consistent results with user-friendly features and robust safety mechanisms.

Model: HX-C83502 HX-C83508

The PowerSpin HX comes with 2 Rotor options:

6 Place Fixed Angle for 2-15ml tubes 8 Place Fixed Angle for 2-15ml tubes






8 Place Fixed Angle

UNICO







182 Ridge Road, Suite E, Dayton, NJ 08810 USA

Email: info@unicosci.com Phone: 732-274-1155



182 Ridge Road, Suite E, Dayton, NJ 08810 USA

Email: info@unicosci.com Phone: 732-274-1155

UNICO® PowerSpin™ LX Series Centrifuge

Maximum Speed: 300-4000 RPM

Maximum RCF: 1860 xg Capacity: 15ml × 8 tubes Noise Level: ≤60 dB(A) Drive Motor: Brushless DC

Power Supply: 100V-240V, 50Hz/60Hz

The UNICO PowerSpin LX Centrifuge combines high performance with versatile functionality, making it an ideal choice for various laboratory applications. Designed with a robust brushless DC motor and user-friendly controls, this centrifuge offers reliable and efficient separation for a range of sample types. It features preset programs for blood, urine, and fecal samples, enhancing convenience and accuracy in laboratory settings.

Model: LX-C85608





Maximum Speed: 3500 rpm Maximum RCF: 1397×g Capacity: 15ml x 8 tubes Noise Level: ≤60 dB(A) Drive Motor: Brushless DC Power Supply: DC 12V, 3A

Elevate your laboratory's efficiency and precision with the UNICO PX Centrifuge, a state-of-the-art instrument designed to meet the demands of diverse scientific and clinical applications. Combining advanced technology with user-friendly features, the UNICO PX delivers exceptional performance.

Model: PX-C82608

UNICO® PowerSpin™ FX Series Centrifuge

Maximum Speed: 3500 rpm
Maximum RCF: 1397×g
Capacity: 15ml x 8 tubes
Noise Level: ≤60 dB(A)
Drive Motor: Brushless DC
Power Supply: 110V, 60Hz, 2A

The FX Centrifuge is a high-performance laboratory instrument designed to deliver accurate and reliable results for a wide range of applications. With its advanced technology and user-friendly features, the FX Centrifuge provides a perfect balance of speed, capacity, and safety. Ideal for clinical, research, and quality control environments, this centrifuge ensures efficient sample processing and optimal results in a compact, durable package.

Model: FX-C80508



UNICO® PowerSpin™ PXV Series Centrifuge

Maximum Speed: 500 to 3800 rpm

Maximum RCF: 1647×g
Capacity: 15ml x 8 tubes
Noise Level: ≤60 dB(A)
Drive Motor: Brushless DC
Power Supply: DC 12V, 3A

Experience cutting-edge centrifugation technology with the UNICO PXV Centrifuge, a high-performance laboratory instrument designed to deliver superior results across a variety of scientific and clinical applications. The PXV combines robust functionality with advanced features, making it an essential tool for precise sample processing in any laboratory setting

Model: PXV-C82708





182 Ridge Road, Suite E, Dayton, NJ 08810 USA

Email: info@unicosci.com Phone: 732-274-1155



182 Ridge Road, Suite E, Dayton, NJ 08810 USA

Email: info@unicosci.com Phone: 732-274-1155

UNICO® VARIABLE SPEED TUBE ROCKERS



Elevate your laboratory's efficiency with the UNICO® L-TTR450 Double Decker Blood Tube Rocker. Designed for heavy-duty performance while maintaining near-silent operation, this rocker ensures gentle mixing of specimens, preserving the integrity of fragile blood cells. Ideal for preparing hematology controls, it guarantees accurate results without compromising sample quality.

Capacity: Available in configurations for 32, 16, or 10 tubes to suit your workflow needs. **Versatile Tube Size Compatibility**: Accommodates a range of tube sizes up to 16x120mm, making it perfect for various applications.

Variable Speed Control: Adjust mixing speed from 10 to 30 RPM (±2 RPM) to achieve optimal results for delicate samples.

Safe and Secure Design: Equipped with a non-slip, reversible, and autoclavable silicon rubber tube pad, ensuring stability during operation.

Item#	Speed	Capacity	Tubes
L-TTR450	10 - 30 RPM	32	2 - 15 ml

UNICO® VARIABLE SPEED TUBE ROCKERS



The UNICO® Variable Speed Tube Rockers are heavy-duty, yet nearly silent rocker designed for years of everyday use. Smooth movement assures a gently mixed specimen while preserving the integrity of fragile blood cells. The Tube Rocker will also quickly get hematology controls ready to use for accurate results. The non-skid platform also accommodates larger, hard to mix tubes.

Item#	Speed	Capacity	Tubes
L-TTR150	10 - 30 RPM	10	2 - 15 ml
L-TTR250	10 - 30 RPM	16	2 - 15 ml