


PowerSpin™ PX Centrifuge


Model: C827

- POWERSPIN™ PORTA-SPIN™ CENTRIFUGE
- VARIABLE SPEED 500-3800 RPM
- 6 PLACE FIXED ANGLE ROTOR
- 30 MIN. MECHANICAL TIMER
- 6X10ML CAPACITY
- POWER INPUT 15VDC



The PortaSpin™ PX allows you to take your practice on the road. Whether you are a mobile vet, visiting nurse or doctor, we have a PortaSpin that will meet your needs. Available with 6 or 8 place rotors and adapters to fit tubes ranging from 2 – 15 ml.

Standard with 12 Volt car lighter adapter, a 110V electrical outlet adapter is optional for use back at the lab. When using this unit you will most likely not be in your lab setting so you will not have your standard bench-top tube rack. We strongly recommend purchasing the optional Tube Holdster, attaches easily to the side of the centrifuge and will hold up to 18 tubes or adapters



ADDITIONAL INFORMATION

ACCESSORIES	Accessories	
	C800-00	8 Place Rotor
	C800-01	8 Place Rotor
	C800-75	Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case
	C800-02	Regular tube shield, pack of 2
	C800-03	Adapter for short tubes, pack of 10
	C800-04	Pediatric or microtube adapter, pack of 2
	C800-05	Tube Lift Adapters, Pk/12
	C800-12	Balance Tube Holdster
	C826-06	Porta-Spin 12 VDC Adapter, Plugs in to Car Charger, Replacement
	C826-12	Porta-Spin Ac Adapter, 110 Volt, USA Power Cord
SPECIFICATION	Number of Tubes	6
	Tube Capacity	2-15 ml
	Blood	Yes
	Specimens Fecal	Not recommended
	Urine	Not recommended
	Microhematocrit	Not recommended
	Speed	Fixed
	Rotor	Fixed Angle
	Maximum RPM	3400
	Maximum RCF(g)	1320
SIZE AND WEIGHT	Size	10.2 in x 11.8 in x 7.9 in
	Shipping Size	13.4 in x 14.8 in x 10.6 in
	Weight	18 lbs
	Shipping Weight	19 lbs

