


PowerSpin™ PX Centrifuge

Model: C828H

- POWERSPIN™ PORTA-SPIN™ CENTRIFUGE
- FIXED SPEED 3700 RPM
- 8 PLACE FIXED ANGLE ROTOR
- 30 MIN. MECHANICAL TIMER
- 8X10ML CAPACITYR WITH 18 PLACE TUBE HOLDSTER™ RACK
- POWER INPUT 15VDC



The PortaSpin™ PX allows you to take your practice on the road. Whether you are a mobile vet, visiting nurse or doctor, we have a PortaSpin that will meet your needs. Available with 6 or 8 place rotors and adapters to fit tubes ranging from 2 - 15 ml.

Standard with 12 Volt car lighter adapter, a 110V electrical outlet adapter is optional for use back at the lab. When using this unit you will most likely not be in your lab setting so you will not have your standard bench-top tube rack. We strongly recommend purchasing the optional Tube Holdster, attaches easily to the side of the centrifuge and will hold up to 18 tubes or adapters



ADDITIONAL INFORMATION

ACCESSORIES	Accessories <ul style="list-style-type: none"> C800-00 8 Place Rotor C800-01 8 Place Rotor C800-75 Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case C800-02 Regular tube shield, pack of 2 C800-03 Adapter for short tubes, pack of 10 C800-04 Pediatric or microtube adapter, pack of 2 C800-05 Tube Lift Adapters, Pk/12 C800-12 Balance Tube Holdster C826-06 Porta-Spin 12 VDC Adapter, Plugs in to Car Charger, Replacement C826-12 Porta-Spin Ac Adapter, 110 Volt, USA Power Cord
SPECIFICATION	<ul style="list-style-type: none"> Number of Tubes 6 Tube Capacity 2-15 ml Blood Yes Specimens Fecal Not recommended Urine Not recommended Microhematocrit Not recommended Speed Fixed Rotor Fixed Angle Maximum RPM 3400 Maximum RCF(g) 1320
SIZE AND WEIGHT	<ul style="list-style-type: none"> Size 10.2 in x 11.8 in x 7.9 in Shipping Size 13.4 in x 14.8 in x 10.6 in Weight 13 lbs Shipping Weight 15 lbs

