

# PowerSpin™ PX Centrifuge

Model: C829H

- POWERSPIN™ PORTA-SPIN™ CENTRIFUGE
- VARIABLE SPEED 500-3800 RPM
- 8 PLACE FIXED ANGLE ROTOR
- 30 MIN. MECHANICAL TIMER
- 8X10ML CAPACITY WITH 18 PLACE TUBE HOLDSTER™ RACK
- POWER INPUT 15VDC



The PortaSpin™ PX allows you to take your practice on the road. Whether you are a mobile vet, visiting nurse or doctor, we have a PortaSpin that will meet your needs. Available with 6 or 8 place rotors and adapters to fit tubes ranging from 2 – 15 ml.

Standard with 12 Volt car lighter adapter, a 110V electrical outlet adapter is optional for use back at the lab. When using this unit you will most likely not be in your lab setting so you will not have you standard bench-top tube rack. We strongly recommend purchasing the optional Tube Holdster, attaches easily to the side of the centrifuge and will hold up to 18 tubes or adapters



## ADDITIONAL INFORMATION

<b>ACCESSORIES</b>	<b>Accessories</b> <ul style="list-style-type: none"> <li>C800-00 8 Place Rotor</li> <li>C800-01 8 Place Rotor</li> <li>C800-75 Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case</li> <li>C800-02 Regular tube shield, pack of 2</li> <li>C800-03 Adapter for short tubes, pack of 10</li> <li>C800-04 Pediatric or microtube adapter, pack of 2</li> <li>C800-05 Tube Lift Adapters, Pk/12</li> <li>C800-12 Balance Tube Holdster</li> <li>C826-06 Porta-Spin 12 VDC Adapter, Plugs in to Car Charger, Replacement</li> <li>C826-12 Porta-Spin Ac Adapter, 110 Volt, USA Power Cord</li> </ul>
<b>SPECIFICATION</b>	<ul style="list-style-type: none"> <li><b>Number of Tubes</b> 8</li> <li><b>Tube Capacity</b> 2-15 ml</li> <li><b>Blood</b> Yes</li> <li><b>Specimens Fecal</b> Not recommended</li> <li><b>Urine</b> Not recommended</li> <li><b>Microhematocrit</b> Not recommended</li> <li><b>Speed</b> Fixed</li> <li><b>Rotor</b> Fixed Angle</li> <li><b>Maximum RPM</b> 3400</li> <li><b>Maximum RCF(g)</b> 1320</li> </ul>
<b>SIZE AND WEIGHT</b>	<ul style="list-style-type: none"> <li><b>Size</b> 10.2 in x 11.8 in x 7.9 in</li> <li><b>Shipping Size</b> 13.4 in x 14.8 in x 10.6 in</li> <li><b>Weight</b> 18 lbs</li> <li><b>Shipping Weight</b> 19 lbs</li> </ul>

