

PowerSpin™ PX Centrifuge

Model: C829H

- POWERSPIN™ PORTA-SPIN™ CENTRIFUGE
- VARIABLE SPEED 500-3800 RPM
- 8 PLACE FIXED ANGLE ROTOR
- 30 MIN. MECHANICAL TIMER
- 8X10ML CAPACITY WITH 18 PLACE TUBE HOLDSTER™ RACK
- POWER INPUT 15VDC



The UNICO C829H PXV 12V Variable Speed Centrifuge is engineered for flexibility and reliability, making it ideal for both laboratory and mobile applications. This advanced centrifuge features an 8-place rotor that efficiently accommodates a range of tube sizes, and comes with an 18-place Tube Holdster Rack for convenient tube storage and handling.

With its variable speed control, the C829H PXV allows precise adjustments from 500 to 3,800 RPM, catering to various centrifugation needs. It operates on a 12V power system, making it perfect for use in vehicles or other environments with limited power sources. The included vehicle power adapter cable ensures seamless operation on the go.

For added safety, the centrifuge is equipped with a dual lid lock system that prevents the lid from being opened during operation. Additionally, a leveling feature maintains stability and ensures accurate performance. Constructed from high-impact plastic, the C829H PXV is built to withstand frequent use and demanding conditions.

Combining variable speed capabilities, comprehensive tube handling with the 18-place rack, and robust safety features, the UNICO C829H PXV 12V Variable Speed Centrifuge offers a dependable and versatile solution for clinical, industrial, and field-based applications.



ADDITIONAL INFORMATION

TOP SELECTION	No
ACCESSORIES	<p>Accessories</p> <p>C800-00 8 Place Rotor</p> <p>C800-01 8 Place Rotor</p> <p>C800-75 Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case</p> <p>C800-02 Regular tube shield, pack of 2</p> <p>C800-03 Adapter for short tubes, pack of 10</p> <p>C800-04 Pediatric or microtube adapter, pack of 2</p> <p>C800-05 Tube Lift Adapters, Pk/12</p> <p>C800-12 Balance Tube Holdster</p> <p>C826-06 Porta-Spin 12 VDC Adapter, Plugs in to Car Charger, Replacement</p> <p>C826-12 Porta-Spin Ac Adapter, 110 Volt, USA Power Cord</p>
SPECIFICATION	<p>Number of Tubes 8</p> <p>Tube Capacity 2-15 ml</p> <p>Blood Yes</p> <p>Specimens Fecal Not recommended</p> <p>Urine Not recommended</p> <p>Microhematocrit Not recommended</p> <p>Speed Fixed</p> <p>Rotor Fixed Angle</p> <p>Maximum RPM 3400</p> <p>Maximum RCF(g) 1320</p>
SIZE AND WEIGHT	<p>Size 10.2 in x 11.8 in x 7.9 in</p> <p>Shipping Size 13.4 in x 14.8 in x 10.6 in</p>



Weight	18 lbs
---------------	--------

Shipping Weight	19 lbs
------------------------	--------

