

# PowerSpin™ HX Centrifuge

Model: C8306



- POWERSPIN™ HXV CENTRIFUGE
- VARIABLE SPEED 500-3500 RPM
- 6 PLACE HORIZONTAL ROTOR
- 99 MIN. DIGITAL TIMER
- 6 X10ML CAPACITY
- POWER INPUT 110V



The UNICO® PowerSpin™ HX Series Centrifuge (Model # C8306) is a highly reliable and versatile laboratory tool, designed to meet the demands of clinical, industrial, and biological labs. It features a 6-place swing-out rotor, perfect for 3-15 ml tubes, ensuring efficient separation of blood and other materials. With easy-to-use push button keypad controls, the centrifuge can be operated with gloves on, making it convenient for busy lab environments.



ADDITIONAL INFORMATION

TOP SELECTION	No
ACCESSORIES	<div>Accessories (Sold Separately)</div> <div><div>C8300-01</div><div>6 place horizontal rotor for HXV C8300 centrifuge</div></div> <div><div>C8300-02</div><div>8 place horizontal rotor for HXV C8300 centrifuge</div></div> <div><div>C8300-03</div><div>12 place angled rotor for HXV C8300 centrifuge</div></div>
SPECIFICATION	<div><div>Speed (RPM)</div><div>500-3500</div></div> <div><div>Maximum RCF (g)</div><div>1777g</div></div> <div><div>Timer</div><div>Digital</div></div> <div><div>Nominal Capacity</div><div>6 x 15 ml</div></div> <div><div>Rotor</div><div>Horizontal</div></div> <div><div>Lock</div><div>ZERO RPM Digital</div></div> <div><div>Motor</div><div>AC</div></div> <div><div>Maximum Test Tube Size</div><div>16 x 110 mm</div></div>
WARRANTY	1 years

