


# C8312 PowerSpin™ HXV Centrifuge

Model: C8312



- POWERSPIN™ HXV CENTRIFUGE
- VARIABLE SPEED 500-3500 RPM
- 12 PLACE FIXED ANGLE ROTOR
- 99 MIN. DIGITAL TIMER
- 12 X10ML CAPACITY
- POWER INPUT 110V



The UNICO® PowerSpin™ HXV Series Centrifuge combines precision, reliability, and advanced safety features to meet the demanding needs of clinical, industrial, and biological laboratories. Designed for long-lasting performance, the HXV series is perfect for applications requiring blood separation and other material processing. The HXV series includes intuitive push-button keypad controls that are easy to use, even while wearing gloves, allowing for quick adjustments to the centrifuge's settings. With a variable speed range of 500 to 3500 RPM, you can fine-tune the speed to suit your specific protocol needs. An integrated electronic timer, adjustable from 1 to 99 minutes, features a beep signal at the end of the run, while the unit's automatic memory reset ensures consistent results across repeated protocols by recalling the last time and speed settings once the centrifuge is powered on again.



## ADDITIONAL INFORMATION

<b>TOP SELECTION</b>	No
<b>ACCESSORIES</b>	<p><b>Accessories (Sold Separately)</b></p> <p>C8300-01                      6 place horizontal rotor for HXV C8300 centrifuge</p> <p>C8300-02                      8 place horizontal rotor for HXV C8300 centrifuge</p> <p>C8300-03                      12 place angled rotor for HXV C8300 centrifuge</p> <p>C800-02                        Regular tube shield, pack of 2</p> <p>C800-03                        Adapter for short tubes, pack of 10</p> <p>C800-04                        Pediatric or microtube adapter, pack of 2</p> <p>C800-05                        Tube Lift Adapters, Pk/12</p> <p>C800-12                        Balance Tube Holdster</p>
<b>SPECIFICATION</b>	<p><b>Speed (RPM)</b>                      500-3500</p> <p><b>Maximum RCF (g)</b>                      1777g</p> <p><b>Timer</b>                                      Digital</p> <p><b>Nominal Capacity</b>                      12 x 15 ml</p> <p><b>Rotor</b>                                      Fixed Angle</p> <p><b>Lock</b>                                        ZERO RPM Digital</p> <p><b>Motor</b>                                        AC</p> <p><b>Maximum Test Tube Size</b>                      16 x 110 mm</p>
<b>WARRANTY</b>	1 years

