

PowerSpin™ LX Centrifuge

Model: C856

- POWERSPIN™ LX CENTRIFUGE
- VARIABLE SPEED 300-4000 RPM
- 6 PLACE FIXED ANGLE ROTOR
- 30 MIN. MECHANICAL TIMER
- 6X10ML OR 3X15ML CAPACITY
- POWER INPUT 110V



UNICO® PowerSpin™ Series centrifuges are compact, economical and quiet. We offer both fixed (FX) and variable (LX) speed models with fixed angle rotors, for 6 or 8 tubes to spin 8- 15 ml tubes. By using adapters, you can also spin tubes ranging from 2 - 7 ml. An optional Tube Holdster™ allows you to have up to 18 tubes or adapters at your fingertips, side-mounted on the centrifuge. PowerSpins are pre-drilled/tapped to accept optional Tube Holdsters.



ADDITIONAL INFORMATION

TOP SELECTION	No
ACCESSORIES	<p>Accessories</p> <p>C800-00 8 Place Rotor</p> <p>C800-01 8 Place Rotor</p> <p>C800-75 Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case</p> <p>C800-02 Regular tube shield, pack of 2</p> <p>C800-03 Adapter for short tubes, pack of 10</p> <p>C800-04 Pediatric or microtube adapter, pack of 2</p> <p>C800-05 Tube Lift Adapters, Pk/12</p> <p>C800-12 Balance Tube Holdster</p>
SPECIFICATION	<p>Number of Tubes 6</p> <p>Tube Capacity 2-15 ml</p> <p>Blood Yes</p> <p>Specimens Fecal Yes</p> <p>Urine Yes</p> <p>Microhematocrit With Adapters</p> <p>Speed Variable</p> <p>Rotor Fixed</p> <p>Maximum RPM 300-4,000</p> <p>Maximum RCF(g) Up to 1,828</p>
SIZE AND WEIGHT	<p>Size 10.2 in x 11.8 in x 7.9 in</p> <p>Shipping Size 13.4 in x 14.8 in x 10.6 in</p> <p>Weight 16 lbs</p> <p>Shipping Weight 18 lbs</p>

