

# PowerSpin™ LX Centrifuge

Model: C858

- POWERSPIN™ LX CENTRIFUGE
- VARIABLE SPEED 300-4000 RPM
- 8 PLACE FIXED ANGLE ROTOR
- 30 MIN. MECHANICAL TIMER
- 6X10ML OR 3X15ML CAPACITY
- POWER INPUT 110V



UNICO® PowerSpin™ Series centrifuges are compact, economical and quiet. We offer both fixed (FX) and variable (LX) speed models with fixed angle rotors, for 6 or 8 tubes to spin 8– 15 ml tubes. By using adapters, you can also spin tubes ranging from 2 – 7 ml. An optional Tube Holdster™ allows you to have up to 18 tubes or adapters at your fingertips, side-mounted on the centrifuge. PowerSpins are pre-drilled/tapped to accept optional Tube Holdsters.



ADDITIONAL INFORMATION

ACCESSORIES	Accessories	
	C800-00	8 Place Rotor
	C800-01	8 Place Rotor
	C800-75	Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case
	C800-02	Regular tube shield, pack of 2
	C800-03	Adapter for short tubes, pack of 10
	C800-04	Pediatric or microtube adapter, pack of 2
	C800-05	Tube Lift Adapters, Pk/12
	C800-12	Balance Tube Holdster
SPECIFICATION	Number of Tubes	8
	Tube Capacity	2-15 ml
	Blood	Yes
	Specimens Fecal	Yes
	Urine	Yes
	Microhematocrit	With Adapters
	Speed	Variable
	Rotor	Fixed
	Maximum RPM	300-4,000
SIZE AND WEIGHT	Maximum RCF(g)	Up to 1,828
	Size	10.2 in x 11.8 in x 7.9 in (260 mm x 300 mm x 200 mm)
	Shipping Size	13.4 in x 14.8 in x 10.6 in (340 mm x 375 mm x 270 mm)
	Weight	16 lbs (7.3 kg)
	Shipping Weight	18 lbs (8.2 kg)

