


# PowerSpin™ LX Centrifuge

Model: C858D

- POWERSPIN™ LX CENTRIFUGE
- DIGITAL VARIABLE SPEED 300 - 4000 RPM
- 8 PLACE FIXED ANGLE ROTOR
- 30 MIN MECHANICAL TIMER
- 6X10ML OR 3X15ML CAPACITY
- POWER INPUT 110V



UNICO® PowerSpin™ Series centrifuges are compact, economical and quiet. We offer both fixed (FX) and variable (LX) speed models with fixed angle rotors, for 6 or 8 tubes to spin 8- 15 ml tubes. By using adapters, you can also spin tubes ranging from 2 - 7 ml. An optional Tube Holdster™ allows you to have up to 18 tubes or adapters at your fingertips, side-mounted on the centrifuge. PowerSpins are pre-drilled/tapped to accept optional Tube Holdsters.

**ADDITIONAL INFORMATION**

<b>ACCESSORIES</b>	<b>Accessories</b>	
	C800-00	8 Place Rotor
	C800-01	8 Place Rotor
	C800-75	Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case
	C800-02	Regular tube shield, pack of 2
	C800-03	Adapter for short tubes, pack of 10
	C800-04	Pediatric or microtube adapter, pack of 2
	C800-05	Tube Lift Adapters, Pk/12
	C800-12	Balance Tube Holdster
<b>SIZE AND WEIGHT</b>	<b>Size</b>	10.2 in x 11.8 in x 7.9 in
	<b>Shipping Size</b>	13.4 in x 14.8 in x 10.6 in
	<b>Weight</b>	16 lbs
	<b>Shipping Weight</b>	18 lbs

