


PowerSpin™ DX Centrifuge

Model: C8706



- POWERSPIN™ DX CENTRIFUGE
- VARIABLE SPEED 500-4000 RPM
- PRE-PROGRAMMED OPTIONS FOR BLOOD, URINE AND FECAL
- 6 PLACE HORIZONTAL ROTOR
- 1:00-99:99 MIN DIGITAL TIMER
- CAPACITY 6X15 ML TUBES
- POWER INPUT 110V



The UNICO PowerSpin™ DX Centrifuge is engineered to meet the demanding needs of busy veterinary clinics. Specifically, the C8706 model, equipped with a 6-place horizontal rotor, is ideal for processing fecal samples, making it perfectly compatible with the VETSCAN IMAGYST™ system for parasite detection and other veterinary diagnostic tests. Featuring digital variable speed control (500-4000 RPM), an electronic timer, and a user-friendly interface with preset programs for common applications, the C8706 ensures efficient and reliable sample processing. Its compact design and quiet operation minimize disruptions in your clinic, while the maintenance-free brushless DC motor ensures long-lasting performance.



NOT INCLUDED



ADDITIONAL INFORMATION

TOP SELECTION	No	
ACCESSORIES	Accessories (Sold Separately) C8600-01 6 Place Horizontal Rotor, 2-15ml for DX Centrifuge C8600-02 4 Place Horizontal rotor, 50ml DX Centrifuge C8600-03M 24 Place Fixed Angle Rotor, 2-15ml for DX Centrifuge, Metal (without Tube Shields)	
SPECIFICATION	Speed (RPM)	500-4000rpm
	Maximum RCF (g)	1848g
	Rotor	6 x 15ml Horizontal
	Tube Size	2-15ml
	Motor	DC
	Power Input	110V
SIZE AND WEIGHT	Size	13.75 inch Height (22 inch height with lid open) x 15 inch Deep x 11 inch Width
	Shipping Size	19"X18"X15"
	Weight	44 LBS
	Shipping Weight	44 LBS
WARRANTY	1 years	

