

# PowerSpin™ DX Centrifuge

Model: C8706E

- POWERSPIN™ DX CENTRIFUGE
- VARIABLE SPEED 1000-3400 RPM
- VARIABLE SPEED 500-3500 RPM
- PRE-PROGRAMMED OPTIONS FOR BLOOD, URINE AND FECAL
- 6 PLACE HORIZONTAL ROTOR
- 60 MIN DIGITAL TIMER
- CAPACITY 6X15 ML TUBES
- POWER INPUT 220V



UNICO is proud to present the PowerSpin™ DX Centrifuge

With 3 rotor options this centrifuge is a great answer for every day laboratory tasks. Choose from 24 place fixed rotor for 2-15ml tubes (#C8724), 6 place horizontal rotor for 2-15ml tubes(#C8706) or 4 place horizontal rotor for up to 50ml tubes (#C8704)

The PowerSpin DX series features a large, spacious stainless steel lined interior, all digital push-button control with 3 preset operational modes for Blood, Urine and Fecal, variable speed 500-4000 RPM and zero RPM lid latching system. DX's are available in three main versions, fixed angle and horizontal spin rotors. All DX's feature maintenance-free brushless motors and suction cup rubber feet.



## ADDITIONAL INFORMATION

<b>ACCESSORIES</b>	<b>Accessories</b> <ul style="list-style-type: none"> <li>C8600-01                      6 Place Horizontal Rotor, 2-15ml for DX Centrifuge</li> <li>C8600-02                      4 Place Horizontal rotor, 50ml DX Centrifuge</li> <li>C8600-03M                    24 Place Fixed Angle Rotor, 2-15ml for DX Centrifuge, Metal (without Tube Shields)</li> <li>C800-02                        Regular tube shield, pack of 2</li> <li>C800-03                        Adapter for short tubes, pack of 10</li> <li>C800-04                        Pediatric or microtube adapter, pack of 2</li> <li>C800-05                        Tube Lift Adapters, Pk/12</li> </ul>
<b>SPECIFICATION</b>	<ul style="list-style-type: none"> <li><b>Speed (RPM)</b>                500~4000 rpm</li> <li><b>Maximum RCF (g)</b>         2508</li> <li><b>Rotor</b>                         6 place horizontal</li> <li><b>Tube Size</b>                    2-15ml</li> <li><b>Motor</b>                         DC</li> <li><b>Power Input</b>                220V</li> </ul>
<b>SIZE AND WEIGHT</b>	<ul style="list-style-type: none"> <li><b>Size</b>                            13.75 inch Height (22 inch height with lid open) x 15 inch Deep x 11 inch Width</li> <li><b>Shipping Size</b>             19"X18"X15"</li> <li><b>Weight</b>                        31 LBS</li> <li><b>Shipping Weight</b>         36 LBS</li> </ul>

