

Basic Lab Microscope

Model: M280

- SEIDENTOPF BINOCULAR HEAD
- WF10X PAIRED EYEPIECES
- 4X,10X,40XR,100XR DIN ACHROMAT OBJECTIVES
- BUILT-IN MECHANICAL STAGE WITH ABBE NA1.25 CONDENSER & IRIS
- LED ILLUMINATION



UNICO M280 Series Microscope is small in size, but big on features. With 4 position reverse nosepiece, mechanical stage great optics and variable LED illumination for years of reliable service.



ADDITIONAL INFORMATION

TOP SELECTION	Yes
ACCESSORIES	Bodies
	M280-3004 SEIDENTOPF BINOCULAR FOR M280 SERIES
	M280-3005 SEIDENTOPF TRINOCULAR FOR M280 SERIES
	Optics
	M280-1401 10X WIDEFIELD EYEPIECE, SINGLE
	M280-2101 4X DIN ACHROMAT OBJECTIVE, N.A. 0.10
M280-2102 10X DIN ACHROMAT OBJECTIVE, N.A. 0.25	



	M280-2103	40X DIN ACHROMAT OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS
	M280-2104	100X DIN ACHROMAT OBJECTIVE, N.A. 1.25, OIL, RETRACTABLE FRONT LENS
	M280-2301	4X DIN PLAN ACHROMAT OBJECTIVE, N.A. 0.10
	M280-2302	10X DIN PLAN ACHROMAT OBJECTIVE, N.A. 0.25
	M280-2303	40X DIN PLAN ACHROMAT OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS
	M280-2304	100X DIN PLAN ACHROMAT OBJECTIVE, N.A. 1.25, OIL, RETRACTABLE FRONT LENS
	Miscellaneous Accessories	
	M280-5003	MECHANICAL STAGE WITH ABBE N.A. 1.25 CONDENSER AND IRIS DIAPHRAGM
	M280-7003	FUSE 250V 1.5A
	M280-8001	DUST COVER FOR M280 SERIES
	M280-8902	EYEPIECE POINTER FOR M280 SERIES
SPECIFICATION	Head	Binocular
	Eyepiece	WF10X
	Objectives	Achromat 4X/10X/40X/100X(oil)
	Condenser	NA 1.25
	Diaphragm	Iris
	Stage	Mechanical
	Illumination	LED
	Bulb	3 Watt LED
	Focusing	Coaxial

