


# PowerSpin™ DX Centrifuge

Model: C8724HE

- POWERSPIN™ DX CENTRIFUGE
- VARIABLE SPEED 1000-3400 RPM
- 24 PLACE FIXED ANGEL ROTOR
- 60 MIN DIGITAL TIMER
- CAPACITY 24X10 ML TUBES OR 16X15 ML TUBES
- 24 PLACE TUBE HOLDSTER™ RACK
- POWER INPUT 110V



UNICO is proud to present the PowerSpin™ DX Centrifuge

With 3 rotor options this centrifuge is a great answer for every day laboratory tasks. Choose from 24 place fixed rotor for 2-15ml tubes (#C8724), 6 place horizontal rotor for 2-15ml tubes(#C8706) or 4 place horizontal rotor for up to 50ml tubes (#C8704)

The PowerSpin DX series features a large, spacious stainless steel lined interior, all digital push-button control with 3 preset operational modes for Blood, Urine and Fecal, variable speed 500-4000 RPM and zero RPM lid latching system. DX's are available in three main versions, fixed angle and horizontal spin rotors. All DX's feature maintenance-free brushless motors and suction cup rubber feet.

**ADDITIONAL INFORMATION**

<b>ACCESSORIES</b>	<b>Accessories</b>
	C8600-01 6 Place Horizontal Rotor, 2-15ml for DX Centrifuge
	C8600-02 4 Place Horizontal rotor, 50ml DX Centrifuge
	C8600-03M 24 Place Fixed Angle Rotor, 2-15ml for DX Centrifuge, Metal (without Tube Shields)
	C800-02 Regular tube shield, pack of 2
	C800-03 Adapter for short tubes, pack of 10
	C800-04 Pediatric or microtube adapter, pack of 2
	C800-05 Tube Lift Adapters, Pk/12
<b>SIZE AND WEIGHT</b>	<b>Size</b> 13.75 inch Height (22 inch height with lid open) x 15 inch Deep x 11 inch Width
	<b>Shipping Size</b> 19"X18"X15"
	<b>Weight</b> 31 LBS
	<b>Shipping Weight</b> 36 LBS

