

PowerSpin™ DX Centrifuge

Model: C8724HE



- POWERSPIN™ DX CENTRIFUGE
- VARIABLE SPEED 1000-3400 RPM
- 24 PLACE FIXED ANGEL ROTOR
- 60 MIN DIGITAL TIMER
- CAPACITY 24X10 ML TUBES OR 16X15 ML TUBES
- 24 PLACE TUBE HOLDSTER™ RACK
- POWER INPUT 110V





ADDITIONAL INFORMATION

TOP SELECTION	No	
ACCESSORIES	Accessories	
	C8600-01	6 Place Horizontal Rotor, 2-15ml for DX Centrifuge
	C8600-02	4 Place Horizontal rotor, 50ml DX Centrifuge
	C8600-03M	24 Place Fixed Angle Rotor, 2-15ml for DX Centrifuge, Metal (without Tube Shields)
	C800-02	Regular tube shield, pack of 2
	C800-03	Adapter for short tubes, pack of 10
	C800-04	Pediatric or microtube adapter, pack of 2
	C800-05	Tube Lift Adapters, Pk/12
SIZE AND WEIGHT	Size	13.75 inch Height (22 inch height with lid open) x 15 inch Deep x 11 inch Width
	Shipping Size	19"X18"X15"
	Weight	31 LBS
	Shipping Weight	36 LBS
WARRANTY	1 years	

