

PowerSpin™ EX Centrifuge

Model: C8924

- POWERSPIN™ EX CENTRIFUGE
- VARIABLE SPEED 500-3500 RPM
- PRE-PROGRAMMED OPTIONS FOR BLOOD, URINE AND FECAL
- 24 PLACE HORIZONTAL ROTOR
- 99 MIN DIGITAL TIMER
- CAPACITY 24X10ML TUBES
- POWER INPUT 110V



UNICO is proud to present the PowerSpin™EX Centrifuge with 24 place horizontal rotor this centrifuge is a great answer for a busy laboratory.

The PowerSpin™ EX has variable speed 500 – 3500 RPM, to insure a clean and complete separation. Quiet and dependable, with a maintenance free motor, with stainless steel tube shields it can be used to process blood, urine, fecal and other specimens.


The all digital front panel with 3 pre-programmed options

3200 RPM for 10 minutes for Blood

1800 RPM for 5 minutes for Urine

1300 RPM for 10 minutes for Fecal

Custom settings are easy and intuitive with all digital panel and oversized display.



ADDITIONAL INFORMATION

ACCESSORIES	Accessories
	C800-03 Adapter for short tubes, pack of 10
	C800-04 Pediatric or microtube adapter, pack of 2
	C800-05 Tube Lift Adapters, Pk/12

