

PowerSpin™ BX Centrifuge

Model: C881



- POWERSPIN™ BX CENTRIFUGE
- VARIABLE SPEED 1000-5000 RPM
- 8 PLACE FIXED ANGLE ROTOR
- 99 MIN DIGITAL TIMER
- CAPACITY 8X10ML OR 8X12ML TUBES
- POWER INPUT 110V



Boost your lab's efficiency with the C881 BX High-Speed Centrifuge. This robust benchtop centrifuge boasts an 8-place rotor, accommodating 8x10ml tubes, and achieves speeds up to 5000 RPM (up to 11,000 RPM with compatible rotors). The oversized LCD display provides clear, real-time data, while the quiet, maintenance-free brushless DC motor ensures reliable operation. A unique feature is the integrated reader on the microhematocrit rotor cover, eliminating the need for separate reading devices. Operates on 110V.



ADDITIONAL INFORMATION

TOP SELECTION	No	
ACCESSORIES	Accessories (Sold Separately) C880-08 Fixed Angle Rotor, 8 place C880-02 Microhematocrit Rotor, 24 place C880-03 Microtube rotor, 24 place C880-10 BX Rotor Rack C881-03 Short tube adapter, for 13x75mm tubes C881-04 Micro / Pediatric tube adapter for 0.5-2.0ml tubes	
SPECIFICATION	Rotor	Fixed Angle, 8 place rotor for 10ml SST Tubes
	Speed Range	1,000 to 5,000 RPM
	Maximum RCF(g)	2571g
	Timer Range	0-99 minutes LCD
	Max. Rotor Capacity	8 x 10ml or 8x12ml
	Dimension	270mm x 320mm x 220mm
SIZE AND WEIGHT	Size	10.5 in x 12.5 in x 9 in (270 mm x 320 mm x 220 mm)
	Shipping Size	19 in x 15 in x 13 in
	Weight	15.5 lbs (C881 & C882); 16.5 lbs (C883)
	Shipping Weight	19 lbs (C881 and C882) ; 20 lbs (C883)
WARRANTY	1 years	

