


# Basic Lab Microscope

Model: M280PL

- SEIDENTOPF BINOCULAR HEAD
- WF10X PAIRED EYEPIECES
- 4X,10X,40XR,100XR DIN PLAN ACHROMAT OBJECTIVES
- BUILT-IN MECHANICAL STAGE WITH ABBE NA1.25 CONDENSER & IRIS
- LED ILLUMINATION



UNICO M280 Series Microscope is small in size, but big on features. With 4 position reverse nosepiece, mechanical stage great optics and variable LED illumination for years of reliable service.



## ADDITIONAL INFORMATION

### ACCESSORIES

#### Bodies

M280-3004	SEIDENTOPF BINOCULAR FOR M280 SERIES
M280-3005	SEIDENTOPF TRINOCULAR FOR M280 SERIES

#### Optics

M280-1401	10X WIDEFIELD EYEPIECE, SINGLE
M280-2101	4X DIN ACHROMAT OBJECTIVE, N.A. 0.10
M280-2102	10X DIN ACHROMAT OBJECTIVE, N.A. 0.25
M280-2103	40X DIN ACHROMAT OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS
M280-2104	100X DIN ACHROMAT OBJECTIVE, N.A. 1.25, OIL, RETRACTABLE FRONT LENS
M280-2301	4X DIN PLAN ACHROMAT OBJECTIVE, N.A. 0.10
M280-2302	10X DIN PLAN ACHROMAT OBJECTIVE, N.A. 0.25
M280-2303	40X DIN PLAN ACHROMAT OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS
M280-2304	100X DIN PLAN ACHROMAT OBJECTIVE, N.A. 1.25, OIL, RETRACTABLE FRONT LENS



	<p><b>Miscellaneous Accessories</b></p> <p>M280-5003                      MECHANICAL STAGE WITH ABBE N.A. 1.25 CONDENSER AND IRIS DIAPHRAGM</p> <p>M280-7003                      FUSE 250V 1.5A</p> <p>M280-8001                      DUST COVER FOR M280 SERIES</p> <p>M280-8902                      EYEPIECE POINTER FOR M280 SERIES</p>	
<p><b>SPECIFICATION</b></p>	<p><b>Head</b>                              Binocular</p> <p><b>Eyepiece</b>                        WF10X</p> <p><b>Objectives</b>                      Plan 4X/10X/40X/100X(oil)</p> <p><b>Condenser</b>                      NA 1.25</p> <p><b>Diaphragm</b>                      Iris</p> <p><b>Stage</b>                              Mechanical</p> <p><b>Illumination</b>                    LED</p> <p><b>Bulb</b>                                3 Watt LED</p> <p><b>Focusing</b>                        Coaxial</p>	

