


PowerSpin™ MH Centrifuge


Model: CMH30



- POWERSPIN™ MH MICRO-HEMATOCRIT CENTRIFUGE
- FIXED SPEED 12000 RPM
- 24 PLACE MICROHEMATOCRIT ROTOR
- 5 MIN DIGITAL TIMER
- POWER INPUT 110V
- MICROHEMATOCRIT READER INCLUDED.



The CMH30 Heavy Duty 24-Place Micro Hematocrit Centrifuge is a reliable and efficient solution for determining hematocrit levels in blood samples. This 110V model is built for demanding laboratory environments and features a robust design capable of handling high volumes. It accommodates up to 24 micro hematocrit capillary tubes, allowing for simultaneous processing and increased throughput. The integrated 5-minute timer ensures consistent and accurate spin cycles, while the convenient braking mechanism allows for quick deceleration and faster turnaround times. This package includes a built-in reader for easy and precise hematocrit determination, eliminating the need for separate equipment. For added convenience and extended use, an extra gasket is also included.



ADDITIONAL INFORMATION

TOP SELECTION	No	
ACCESSORIES	Accessories	
	CMH30-Reader	Microhematocrit reader
	CMH30-03	Gasket Replacement For CMH30 Centrifuge
SPECIFICATION	Speed	12000 RPM
	RCF	14,837g
	Rotor	24 place
	Power Input	110V
SIZE AND WEIGHT	Size	9in x 10 in x 9 in
	Shipping Size	13 in x 13 in x 13 in
	Weight	18 lbs
	Shipping Weight	19.4 lbs
WARRANTY	1 years	

