


Tube Rocker

Model: L-TTR200E

- UNICO TEST TUBE ROCKER
- UP TO 15 TUBES/VIALS CAPACITY
- EASY OPERATION TOP PAUSE BUTTON
- FIXED SPEED OF 25 ROCKING CYCLES PER MINUTE
- POWER INPUT 220 V



The UNICO L-TTR200E Blood Tube Mixer is a robust and efficient solution designed for gentle and effective mixing of blood tubes and vials. With the capacity to handle up to 15 tubes or vials simultaneously, this mixer is ideal for busy laboratories requiring reliable performance for sample preparation. Operating at a consistent 25 rocking cycles per minute, the L-TTR200E ensures thorough mixing while preserving the integrity of sensitive specimens. The mixer features a smooth rocking motion, rocking up to 30 degrees in each direction, which helps in maintaining the quality of blood cells and other delicate samples. Designed for 220V power supplies, the L-TTR200E is built with durability in mind. Its non-slip, reversible silicon rubber pad is autoclavable for easy cleaning and hygiene. Compact and space-efficient, this mixer fits comfortably on benchtops, making it a practical addition to any lab. Whether you are handling blood samples, hematology controls, or other sensitive vials, the UNICO L-TTR200E provides a dependable and gentle mixing solution, enhancing laboratory efficiency and reliability.



ADDITIONAL INFORMATION

| | | |
|------------------------|------------------------|---------------------------|
| TOP SELECTION | No | |
| SPECIFICATION | Capacity | 15 tubes |
| | Speed | 25cpm |
| | Voltage | 220V |
| SIZE AND WEIGHT | Size | 8.6 in x 4.7 in x 5.0 in |
| | Shipping Size | 10.6 in x 6.7 in x 6.9 in |
| | Weight | 4.0 lbs (1.8 kg) |
| | Shipping Weight | 5.0 lbs (2.3 kg) |

