


Vortex Mixer

Model: L-VM1000



- UNICO VORTEX MIXER
- SINGLE TUBE HEAD
- 3-POSITION SWITCH FOR CONSTANT OR "TOUCH-ON" OPERATION
- VARIABLE SPEED UP TO 3000RPM
- POWER INPUT 120V



UNICO economical vortex mixers are rugged and sturdy for your most demanding applications. Dual mode switch lets you vortex continuously or intermittently using the touch mode - simply press your tube to the soft rubber mixing cup. Both model feature speed control to adjust the vortexing velocity to your requirements.



ADDITIONAL INFORMATION

TOP SELECTION	No
SIZE AND WEIGHT	Size 10.4 in x 9.2 in x 6.0 in (263 x 234 x 153 mm) Shipping Size 12.2 in x 11.2 in x 7.4 in (311 x 285 x 187 mm) Weight 9.8 lbs (4.4 kg) Shipping Weight 10.8 lbs (4.9 kg)
WARRANTY	1 years

