

Vortex Mixer

Model: LVM2000E

- UNICO VORTEX MIXER
- 3.75 INCHES VIBRATION PAD
- 3-POSITION SWITCH FOR CONSTANT OR "TOUCH-ON" OPERATION
- VARIABLE SPEED UP TO 3000RPM
- POWER INPUT 220V



UNICO economical vortex mixers are rugged and sturdy for your most demanding applications. Dual mode switch lets you vortex continuously or intermittently using the touch mode - simply press your tube to the soft rubber mixing cup. Both model feature speed control to adjust the vortexing velocity to your requirements.

ADDITIONAL INFORMATION

SIZE AND WEIGHT	Size	10.4 in x 9.2 in x 6.0 in (263 x 234 x 153 mm)
	Shipping Size	12.2 in x 11.2 in x 7.4 in (311 x 285 x 187 mm)
	Weight	9.8 lbs (4.4 kg)
	Shipping Weight	10.8 lbs (4.9 kg)

