

Vortex Mixer

Model: LVM2000E

- UNICO VORTEX MIXER
- 3.75 INCHES VIBRATION PAD
- 3-POSITION SWITCH FOR CONSTANT OR "TOUCH-ON" OPERATION
- VARIABLE SPEED UP TO 3000RPM
- POWER INPUT 220V



UNICO economical vortex mixers are rugged and sturdy for your most demanding applications. Dual mode switch lets you vortex continuously or intermittently using the touch mode - simply press your tube to the soft rubber mixing cup. Both model feature speed control to adjust the vortexing velocity to your requirements.

ADDITIONAL INFORMATION**SIZE AND WEIGHT**

| | |
|------------------------|---|
| Size | 10.4 in x 9.2 in x 6.0 in (263 x 234 x 153 mm) |
| Shipping Size | 12.2 in x 11.2 in x 7.4 in (311 x 285 x 187 mm) |
| Weight | 9.8 lbs (4.4 kg) |
| Shipping Weight | 10.8 lbs (4.9 kg) |

