


Orbital Shaker

Model: L-RT30C



UNICO L-RT30C Digital Rotator is a compact, durable, variable speed orbital mixer that will accommodate a wide variety of printed well slides, culture plates, flasks and small test tube racks, as well as RPR and VDRL slides or 96 well microtiter (immunoassay) plates. .



ADDITIONAL INFORMATION

SIZE AND WEIGHT	Size	10.4 in x 9.2 in x 6.0 in (263 x 234 x 153 mm)
	Shipping Size	12.2 in x 11.2 in x 7.4 in (311 x 285 x 187 mm)
	Weight	9.8 lbs (4.4 kg)
	Shipping Weight	10.8 lbs (4.9 kg)

