

Compact Incubator

Model: L-CU60



- UNICO COMPACT INCUBATOR
- 6L CHAMBER CAPACITY
- REMOVABLE SHELF
- TEMPERATURE AMBIENT TO 60°C (140°F)
- POWER INPUT 110 VOLT



The L-CU60 is a highly accurate incubator that is small in size but big in performance. Featuring a see-through door that makes a visual check on your cultures quick and easy. With a capacity of 10 standard size culture plates per shelf, the nominal capacity is 20 plates. The acrylic body is a great insulator and is easy to clean and maintain. With a small footprint, this unit is perfect for labs with limited benchtop space



ADDITIONAL INFORMATION

TOP SELECTION	No	
SPECIFICATION	Chamber Capacity	6 Liters (0.21 cubic feet)
	Chamber Size: Internal Dimensions	6.8
SIZE AND WEIGHT	Shipping Size	14 in x 13 in x 13 in
	Weight	10 lbs
	Shipping Weight	11 lbs
WARRANTY	1 years	

