


Clinical Incubator

Model: L-CU100

- UNICO CLINICAL INCUBATOR
- 10L CHAMBER CAPACITY
- 2 REMOVABLE SHALVES
- SINGLE SEE-THROUGH DOOR DESIGN
- TEMPERATURE AMBIENT TO 70° C
- POWER INPUT 110 VOLT



A compact incubator, perfect for limited bench top space. Capable of heating from ambient up to 70°C (155°F) quickly and holding a constant temperature. Featuring a single door with a window, so you can visually inspect your cultures.



ADDITIONAL INFORMATION

SPECIFICATION

Chamber Capacity	10 Liters (0.35 cubic feet)
Chamber Size	Internal Dimensions 11.2"W x 7.8"D x 8.6"H
Calibration	Digital control & thermometer access port
Temperature Range	Ambient to 70°C (Ambient to 155° F)
Shelves	Two stainless steel shelves
Temperature Control	Digital Control
Display	Digital Display
Power/Consumption	1.4 Amps / 150 Watts
Heating	Convection
Stability	Will hold constant temperature +/- 0.5° C

