

# Clinical Laboratory Incubator

Model: L-CU200

- UNICO CLINICAL LABORATORY INCUBATOR
- 20L CHAMBER CAPACITY
- 2 REMOVABLE SHALVES
- SINGLE SEE-THROUGH DOOR DESIGN
- TEMPERATURE AMBIENT TO 70° C
- POWER INPUT 110 VOLT



UNICO incubators are designed for plates, tubes (slants), urine and general microbiology applications. The spacious stainless steel-lined interiors can also be used as fluid warmers. These incubators feature a double-door design, solid exterior door and a see-thru inner acrylic door.



## ADDITIONAL INFORMATION

### SPECIFICATION

<b>Chamber Capacity</b>	20 Liters (0.7 cubic feet)
<b>Chamber Size</b>	Internal Dimensions 12.3"W x 12.1"D x 9"H
<b>Calibration</b>	Digital control & thermometer access port
<b>Temperature Range</b>	Ambient to 70°C (Ambient to 155° F)
<b>Shelves</b>	Two removable shelves, perforated for air flow
<b>Temperature Control</b>	Digital Control
<b>Display</b>	Digital Display
<b>Power/Consumption</b>	1.4 Amps / 150 Watts
<b>Heating</b>	Convection
<b>Stability</b>	Will hold constant temperature +/- 0.5° C
<b>Power input</b>	Available in 110V and 220V

