


Clinical Laboratory Incubator

Model: L-CU300E

- UNICO CLINICAL LABORATORY INCUBATOR
- 30L CHAMBER CAPACITY
- 2 REMOVABLE SHALVES
- DOUBLE DOOR DESIGN
- TEMPERATURE AMBIENT TO 70° C
- POWER INPUT 220 VOLT



UNICO incubators are designed for plates, tubes (slants), urine and general microbiology applications. The spacious stainless steel-lined interiors can also be used as fluid warmers. These incubators feature a double-door design, solid exterior door and a see-thru inner acrylic door.



ADDITIONAL INFORMATION

| | | |
|---------------|---------------------|--|
| SPECIFICATION | Chamber Capacity | 30 Liters (1.05 cubic feet) |
| | Chamber Size | Internal Dimensions 13"W x 12.4"D x 12.6"H |
| | Calibration | Digital control & thermometer access port |
| | Temperature Range | Ambient to 70°C (Ambient to 155° F) |
| | Shelves | Two removable shelves, perforated for air flow |
| | Temperature Control | Digital Control |
| | Display | Digital Display |
| | Power/Consumption | 1.4 Amps / 150 Watts |
| | Heating | Convection |
| | Stability | Will hold constant temperature +/- 0.5° C |
| | Power Input | Available in 110V and 220V |

