


Clinical Laboratory Incubator

Model: L-CU300E

- UNICO CLINICAL LABORATORY INCUBATOR
- 30L CHAMBER CAPACITY
- 2 REMOVABLE SHALVES
- DOUBLE DOOR DESIGN
- TEMPERATURE AMBIENT TO 70° C
- POWER INPUT 220 VOLT



UNICO incubators are designed for plates, tubes (slants), urine and general microbiology applications. The spacious stainless steel-lined interiors can also be used as fluid warmers. These incubators feature a double-door design, solid exterior door and a see-thru inner acrylic door.



ADDITIONAL INFORMATION

SPECIFICATION

Chamber Capacity	30 Liters (1.05 cubic feet)
Chamber Size	Internal Dimensions 13"W x 12.4"D x 12.6"H
Calibration	Digital control & thermometer access port
Temperature Range	Ambient to 70°C (Ambient to 155° F)
Shelves	Two removable shelves, perforated for air flow
Temperature Control	Digital Control
Display	Digital Display
Power/Consumption	1.4 Amps / 150 Watts
Heating	Convection
Stability	Will hold constant temperature +/- 0.5° C
Power Input	Available in 110V and 220V

