


Productivity Series Spectrophotometer


Model: S-1200



- UNICO 1200 SERIES VISIBLE SPECTROPHOTOMETER
- 325-1000 NM WAVELENGTH RANGE
- 5NM SPECTRAL BANDWIDTH
- SAMPLE HOLDERS FOR 10 MM SINGLE SQUARE CUVETTE AND V-TYPE HOLDER FOR ROUND TUBE
- USB PORT FOR OPTIONAL SOFTWARE
- INCLUDED ARE 2PCS OF 10MM SQUARE GLASS CUVETTES; 12 PCS OF 10 MM TEST TUBES; DUST COVER; USER MANUAL
- VOLTAGE PRESET AT 110V



The MODEL 1200 Multi-Task Visible Spectrophotometer offers versatile performance for a range of applications. This unit features a 5nm bandpass, accommodating various analytical needs. It includes both a round tube cell holder and a single square cuvette holder for flexible sample handling. Also included is a box of 12 test tubes with a 100mm pathlength. This model is preset at 110V.



ADDITIONAL INFORMATION

TOP SELECTION	No
ACCESSORIES	<div><div><div>Software</div><div>S-1200-401</div><div>Unico S-1200 Spectrophotometer Software. Win 7 Or Win 10 Required. Functions include Standard Curve, Abs./%T/Conc. and Abs. vs. Time. includes USB Cable.</div></div><div><div>S-1205-401</div><div>Unico S-1205 Spectrophotometer Software. Win 7 Or Win 10 Required. Functions include Standard Curve, Abs./%T/Conc., and Kinetics and Scanning. includes USB Cable.</div></div></div> <div><div>Accessories</div><div>S-1200-101</div><div>Test Tube V-Type Cell Holder with Base for 8~25 Mm Diameter Test Tubes for S-1200 Series</div><div>S-1200-102</div><div>Single Square Cuvette Holder with Base for 10 Mm Cuvettes for S-1200 Series</div><div>S-1200-103</div><div>4-Position Cell Holder for 10 Mm Square Cuvettes (Standard in Model 1201) for S-1200 Series</div><div>S-1200-104</div><div>Long Path 4-Cell Holder for Up To 50 Mm Rectangular Cuvettes for S-1200 Series</div><div>S-1200-105</div><div>Water-Jacketed Cell Holder W/Base and Panel for Temperature Controlled Application for S-1200 Series</div><div>S-1200-106</div><div>Universal Base for V-Type Holder, Single Cell Holder and Water-Jacketed Holder for S-1200 Series</div><div>S-1200-111</div><div>Black Block To Set 0%T</div><div>S-1200-512</div><div>Dust Cover for S-1200 Spectrophotometer</div><div>S-1205-515</div><div>Halogen Bulb with Vertical Filament G4 12v20w for S-1205</div><div>S-1200-505</div><div>Halogen Bulb, Package Of 1 (6v/10w, G4 Pin)</div></div>



	<div>Wavelength Accuracy</div> <div>± 2 nm</div>
SIZE AND WEIGHT	<div>Size</div> <div>16 x 12 x 8 inch</div> <div>(406 x 305 x 203mm)</div>
	<div>Shipping Size</div> <div>19.5 x 16.1 x 11.3 inch</div> <div>(496 x 410 x 286 mm)</div>
	<div>Weight</div> <div>18 lbs (8.2 kg)</div>
	<div>Shipping Weight</div> <div>21 lbs (9.5 kg)</div>
WARRANTY	<div>1 years</div>

